ALBUMINA

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Doumas, B.T., Watson, W.A., Biggs, H.G. (1971). Clin. Chim. Acta, 31, 87-96.

ALBUMIN

BROMOCRESOL GREEN METHOD
For "in vitro" determination of Albumin in serum or plasma

Principle

At a certain pH value, albumin is specifically combined with bromocresol green, to produce a coloured complex, which is photometrically measured.

Reagents

Kit 2 x 250 ml. (Ref. 99 72 83). Contents:

A. 2 x 250 ml Bromocresol green. Ref. 99 01 62
B. 1 x 5 ml Standard. Ref. 99 02 46

Kit 3 x 100 ml. (Ref. 99 72 58). Contents:

A. 3 x 100 ml Bromocresol green. Ref. 99 94 82
B. 1 x 5 ml Standard. Ref. 99 02 46

Concentration in the reagent solution are:

Succinate buffer pH 4.2 50 mM
Bromocresol green 0.75 g/L
Surfactants
Preservatives and stabilizers

Standard: Aqueous solution, equivalent to 5 g/dl (50 g/L). Ready to use.

Storage and stability

When stored at 2 -8° C, the reagent will remain stable until the expiration date stated on the label.

Sample

Serum or plasma. Sample may be stored for 2 weeks at 2-8°C, or 4 months at -20 °C.

Caution

The reagent contains Sodium azide at 0.09%. Handle with care. The disposal of the residues has to be made according to legal local regulations.

Performance Characteristics

Assay is linear up to 6 g of Albumin/dl.
The analytical performance characteristics of the product depend both of the reagent and the reading system used, manual or automatic. The following data have been obtained manually.

Intraseres Variation Coefficient: 0.77%
Interseries Variation Coefficient: 0.98%
Recovery: 99.1 %.

Colour stability

While the colour change upon binding of albumin to bromocresol green is almost instantaneous, the colour change from binding with other protein fractions occurs continuously with time, so it is recommended not to delay the absorbance reading.

Autoanalyzers

Technical bulletins for different analyzers are available upon request.

Quality control

Seriscann Normal (Normal Control Serum) (Ref. 99 41 48) and Seriscann Anormal (Abnormal Control Serum) (Ref. 99 46 85).

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Mix well and let stand for 5 min. at room temperature (20-25°C). Read immediately.

Reading
Wavelength: 630 nm.
Blank: BL contents

Calculation
$$\frac{\text{SA Abs.}}{\text{ST Abs.}} \times 5 = \text{g de albumin/dl}$$

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