AMILASA LIQUIDA

METODO IFCC

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
For "in vitro" determination of α -Amylase

Principle

* 4,6-Ethylidene(G_7)-p-nitrophenyl(G_1)- α ,D maltoheptaoside.
* PNP = p-Nitrophenol

Reagents

Kit 1 x 60 ml. (Ref. 99 60 66). Contents:

A. 1 x 50 ml. Buffer / Enzyme. Ref. 99 60 88
B. 1 x 10 ml. Buffer / Substrate. Ref. 99 60 90

Kit 1 x 120 ml. (Ref. 99 26 40). Contents:

A. 1 x 100 ml. Buffer / Enzyme. Ref. 99 26 45
B. 1 x 20 ml. Buffer / Substrate. Ref. 99 26 50

Working reagent

Reagents **A** and **B** are ready-to-use. If a **Monoreagent** procedure is preferred, then the reagents must be mixed in the ratio: 5 parts of **A** (Buffer/Enzyme) + 1 part of **B** (Buffer / substrate).

The concentrations in the reagent solution are:

Hepes buffer pH 7.1 65 mM
NaCl 95 mM
MgCl₂ 15 mM
Ethylidene-G₇PNP 17 mM
 α -Glucosidase ≥ 4500 U/L

Stabilizers and preservatives

Storage and stability

When stored at 2-8°C the components of the kit will remain stable until the expiration date stated on the label. The Monoreagent is stable for 2 weeks at 2-8°C and for 5 days at room temperature ($\leq 25^\circ\text{C}$), when protected from the sunlight.

Sample

Serum or plasma with heparin or urine. Samples free from hemolysis should be used.

In absence of bacterial contamination the α -Amylase activity in serum or plasma will remain stable for 6 days when stored at 2 - 8° C. α -Amylase activity in urine will remain stable for 10 days at 2-8°C, and for 2 days at room.

Procedure	Monoreagent	Bireagent
1.Temp. (37°C.)	sample	sample
	ml	ml
Monoreag. (A+B)	3.00	---
Reagent A	---	2.50
Sample	0.10	0.10
Reagent B	---	0.50

Mix and incubate for 3 min., then read the absorbance and start the stopwatch at the same time. Read the Abs. again, after exactly 1, 2 and 3 min.

Reading

Wavelength: Hg 405 nm. (405 - 420 nm).

Blank: against air.

Cuvette: Thermostated 1 cm light-path.

Calculations

Determine the $\Delta E/\text{min.}$ for every reading and find the mean value.

Serum, plasma, urine

U/L (37° C) = 2930 x ΔE 405 nm./min.

$\mu\text{kal/L}$ (37° C) = 48.8 x ΔE 405 nm./min.

Normal values

Serum / plasma: 28 - 100 U/L (0,47- 1.67 $\mu\text{kal/L}$)

Urine: ≤ 460 U/L (≤ 7.67 $\mu\text{kal/L}$)

Caution

The reagent contains Sodium azide at 0.09%. Handle with care. The disposal of the residues has to be made according to legal local regulations.

Semimicromethod

It can be carried out with half of the volumes as per the macro method.

Performance Characteristics

Linearity: Up to $\Delta E/\text{min.}$ = 0.300. For higher values, it is recommended to dilute the sample 1/10 in saline (NaCl 0.9%) and assay once again.

Multiply the final result by 10.

The analytical performance characteristics of the product depend both of the reagent and the reading system used, manual or automatic. The following data have been obtained manually.

Intraserie Variation Coefficient: 1.56%

Interserie Variation Coefficient: 1.98%

Recovery: 98.8 %.

Sweat and saliva contamination in the glassware will alter the final results.

A slight yellowing of the reagent solution does not affect the determination.

Bilirubin concentrations higher than 15 mg/dl, and Haemoglobin higher than 450 mg/dl will interfere with the reaction.

Quality control

Seriscann Normal (Normal Control Serum) (Ref. 99 41 48) and Seriscann Abnormal (Abnormal Control Serum) (Ref. 99 46 85).

Autoanalyzers

Technical bulletins for different analyzers are available upon request.

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