# CALCIO-ARSENATO III

## METODO COLORIMETRICO

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# CALCIUM-ARSENATO III

## COLORIMETRIC METHOD

For "in vitro" determination of Calcium in serum, plasma or urine

### Principle

At pH 8.5 **Arsenato III** (2,2'-[1,8-Dihydroxy-3,6-disulpho-2,7-naphthalene-bis-(azo)] dibenzearsonic acid) reacts with calcium to form a coloured complex. The intensity developed is proportional to the calcium concentration in the sample.

### Reagents

**Kit 2 x 100 ml. (Ref. 99 24 80).** Contents:

A. 2 x 100 ml. Arsenato III reagent.

Ref. 99 86 10

B. 1 x 5 ml Standard.

Ref. 99 87 19

Aqueous solution of Calcium equivalent to 10 mg/dl (2.49 mmol/L).

Ready-to-use.

The reagent's composition is as follows:

Borate buffer pH 8.5	50 mM
Arsenato III	0.1 mM
8-Hydroxyquinolein	12 mM
Tensoactivos	
Preservatives and stabilizers	

### Storage and stability

The components of the kit, when stored at room temperature (15 - 25°C), will remain stable until the expiration date stated on the label. The reagent and standard are ready-to-use.

### Sample

Serum or plasma with heparin or urine. Serum calcium remains unchanged for 10 days at room temperature (15 -25°C) or at 2-8°C.

When testing a urine sample, it should be previously diluted with deionized water (1 + 1).

Multiply the final result by 2.

### Caution

Handle with care. The disposal of the residues has to be made according to legal local regulations.

Procedure	BL	ST	SA
	ml	ml	ml
Sample	--	--	0.01
Standard	--	0.01	--
Reagent	1.00	1.00	1.00

Mix well and let stand for 3 min. at room temperature (20-25°C).

**Reading**  
Wavelength: 650 nm.  
Blank: BL contents.  
Colour stability: a minimum of 1 hour.

**Calculations**  
SA Abs.  
\_\_\_\_\_ x 10 = mg of calcium / dl  
ST Abs.

**SI Units**  
(mg/dl) x 0.2495 = mmol/L

**Normal values**  
Newborn: 8.0 - 13.0 mg/dl  
Children: 8.8 - 12.0 mg/dl  
Adults: 8.8 - 10.5 mg/dl  
Urine: 100 - 300 mg/24h.

### Performance Characteristics

Reaction is linear up to 15 mg/dl (3.75 mmol/L). If higher concentrations are to be run dilute the sample with deionized water (1 + 1). Multiply the final result by 2. The analytical performance characteristics of the product depend both of the reagent and the reading system used, manual or automatic. The following data have been obtained manually.

Intraserie Variation Coefficient: 2.51%  
Interseries Variation Coefficient: 2.83%  
Recovery: 96.3 %.

Bilirubin concentrations higher than 20 mg/dl and phosphates higher than 40 mg/dl, will interfere with the assay. Those anticoagulants which can behave as chelating agents shall never be used with this test (EDTA, fluoride, oxalate...).

To avoid undesirable contaminations it is expressly recommended the use of disposal plasticware.

### Autoanalyzers

Technical adaptation for different autoanalyzers are available upon request.

### Quality control

Seriscann Normal (Normal Control Serum) (Ref. 99 41 48) and Seriscann Anormal (Abnormal Control Serum) (Ref. 99 46 85).

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