COLESTEROL LÍQUIDO

METODO CHOD – POD

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Cálculos
Abs. PR
_____ x 200 = mg Colesterol/dl
Abs. ST

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
CHOLESTEROL LIQUID

CHOD - POD METHOD

For “in vitro” determination of Cholesterol in serum or plasma.

Principle

Reagents

Kit 1 x 100 ml. (Ref. 99 52 82). Contents:

A. 1 x 100 ml. Reagent. Ref. 99 52 20
B. 1 x 5 ml. Standard. Ref. 99 02 48

Kit 3 x 100 ml. (Ref. 99 52 80). Contents:

A. 3 x 100 ml. Reagent. Ref. 99 52 20
B. 1 x 5 ml. Standard. Ref. 99 02 48

Kit 2 x 250 ml. (Ref. 99 50 12). Contents:

A. 2 x 250 ml. Reagent. Ref. 99 01 59
B. 1 x 5 ml. Standard. Ref. 99 02 48

Working reagent

The reagent is ready-to-use.

Concentrations in the reagent solution are:

Mes buffer pH 6.5 75 mM
Phenol 6 mM
3,5-Dichlorophenol 0,2 mM
4-Aminoantipyrine 0,5 mM
Cholesterol Esterase ≥ 500 kU/L
Cholesterol Oxidase ≥ 300 kU/L
Peroxidase ≥ 1200 kU/L
Non reactive Stabilizers

Standard. Solution of Cholesterol in isopropanol/water equivalent to 200 mg/dl (5,18 mmol/L).

Storage and stability

The components of the kit, stored at 2-8°C, will remain stable until the expiration date stated on the label.

Sample

Serum, or plasma. Samples are stable 1 week at 2-8°C, and up to three months at -20° C.

Procedure	BL	SA	ST
	ml	ml	ml
Sample	--	0.01	--
Standard	--	--	0.01
Working reagent	1.00	1.00	1.00

Mix well and let stand for 5 min. at 37°C, or 10 min at room temperature.

Reading
Wavelength: Hg 546 nm; 505 nm.
Blank: the contents of BL.
Colour stability: 1 hour. (when protected from direct sunlight).

Calculations
SA Abs.
_____ x 200 = mg Cholesterol/dl
ST Abs.

SI Units
(mg/100 dl) x 0.0259 = mmol/L

Clinical interpretation

According to the recommendations of the European Atherosclerosis Society:

Lipid disorder	
Colesterol < 200 mg/dl Triglycerides < 200 mg/dl	NO
Colesterol 200 - 300 mg/dl Colesterol-HDL < 35mg/dl	SI
Colesterol > 300 mg/dl Triglycerides > 200 mg/dl	SI

Caution

The reagent contains phenol. Handle with care. The disposal of the residues has to be made according to legal local regulations.
Do not pipette directly from the bottle of reagent, to avoid undesirable contamination.

Performance Characteristics

The assay is linear up to 700 mg of Cholesterol/dl. For higher concentrations dilute the sample 1/2 with saline (NaCl 0,9%). Multiply the final result by 2.
The analytical performance characteristics of the product depend both of the reagent and the reading system used, manual or automatic. The following data have been obtained manually.

Intraseres Variation Coefficient: 0,87%
Interseries Variation Coefficient: 1,44%
Recovery: 98,6 %.

Hemoglobin concentrations higher than 200 mg/dl interfere the assay, as well as Bilirubin concentrations higher than 15 mg/dl.

Quality control

Seriscann Normal (Normal Control Serum) (Ref. 99 41 48) and Seriscann Abnormal (Abnormal Control Serum) (Ref. 99 46 85).

Autoanalyzers

Technical bulletins for different analyzers, are available upon request.

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