CREATININA

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CREATININE

MODIFIED JAFFE METHOD

For "in vitro" determination of Creatinine in serum, plasma or urine

Principle

At alkaline pH values, creatinine reacts with picric acid to produce a coloured compound, creatinine alkaline picrate, which is photometrically measured.

Reagents

Kit 2 x 100 ml. (Ref. 99 88 91). Contents:

A. 1 x 100 ml. Picric acid solution	Ref. 99 90 06
B. 1 x 100 ml. Alkaline solution	Ref. 99 56 35
C. 1 x 5 ml. Standard.	Ref. 99 83 99

Kit 4 x 250 ml. (Ref. 99 00 77). Contents:

A. 2 x 250 ml. Picric acid solution	Ref. 99 00 82
B. 2 x 250 ml. Alkaline solution	Ref. 99 01 08
C. 1 x 5 ml. Standard.	Ref. 99 83 99

Aqueous solution equivalent to 2 mg/dl, (176.8 µmol/L).
Ready-to-use.

Working reagent

Mix equal parts of each reagent (A and B) prior to assay.
The concentrations in the working reagent are:

Picric acid	55 mM
Sodium carbonate	50 mM
NaOH	0.40 M

Preservatives and stabilizers

Storage and stability

The components of the kit, when stored at room temperature (≤ 25°C), will remain stable until the expiration date stated on the label.

Working reagent is stable for up to 15 days at room temperature (≤ 25°C).
Store in the dark.

Sample

Serum or plasma with heparin and urine. Serum or plasma samples are stable for 24 hours at 2 - 8° C.
To assay urinary creatinine, dilute the sample 1/20 in deionized water.
Multiply the final result by 20.

Caution

The reagent contains Sodium azide at 0.09%. Handle with care. The disposal of the residues has to be made according to legal local regulations.

To calculate the Clearance:

$$\frac{(\text{mg creatinine/dl urine}) \times \text{ml urine}/24 \text{ h.}}{(\text{mg creatinine/dl serum}) \times 1440} = \text{ml/min.}$$

S.I. Units

(mg/dl) x 88.4 = µmol/L

Normal values

	Serum, plasma	Urine
Men:	0.6-1.1 mg/dl	21-26 mg/Kg/24 h.
Women:	0.5-0.9 mg/dl	16-22 mg/Kg/24 h.

Clearance

Men:	97 - 137 ml/min.
Women:	88 - 128 ml/min.

Performance Characteristics

Linearity: Up to 15 mg/dl of Creatinine. For higher values, dilute the sample 1/2 in deionized water and assay once again. Multiply the final result by 2.
The analytical performance characteristics of the product depend both of the reagent and the reading system used, manual or automatic. The following data have been obtained manually.

Intraseres Variation Coefficient: 1.72%
Interseries Variation Coefficient: 2.11%
Recovery: 97.4 %.

Lipemic and hemolyzed sera will interfere with the assay.
Several medicinal substances, eventually present in the sample, such as ascorbic acid or levodopa, can interfere with the final result.
At the described assay conditions, protein and glucose interferences are minimized.

Quality control

Seriscann Normal (Normal Control Serum) (Ref. 99 41 48) and Seriscann Anormal (Abnormal Control Serum) (Ref. 99 46 85).

Autoanalyzers

Technical bulletins for different analyzers are available upon request.

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Procedure	BL	SA	ST
	ml	ml	ml
Standard	--	--	0,1
Sample	--	0,1	--
Working Reagent	1,0	1,0	1,0

Mix and start the stopwatch. Transfer to the measuring cuvette. Read the absorbances at 20 and 80 seconds.

Reading

Wavelength: Hg 546 nm; 510 nm.
Blank: the contents of BL.

Calculations

Determine the ΔAbs. for sample and standard.

$$\Delta \text{Abs.} = \text{Abs.}_{80 \text{ sec.}} - \text{Abs.}_{20 \text{ sec.}}$$
$$\frac{\Delta \text{Abs. SA}}{\Delta \text{Abs. ST}} \times 2 = \text{mg creatinine / dl}$$

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