# FOSFATASA ALCALINA LIQUIDA

## METODO IFCC

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
# ALKALINE PHOSPHATASE LIQUID

## IFCC METHOD

For "in vitro" determination of Alkaline Phosphatase in serum or plasma

### Principle

### Reagents

Kit 1 x 50 ml. (Ref. 99 62 65). Contents:

A. 1 x 40 ml. Buffer solution. Ref. 99 71 04  
B. 1 x 10 ml. Substrate. Ref. 99 84 39

Kit 1 x 125 ml. (Ref. 99 55 15). Contents:

A. 1 x 100 ml. Buffer solution. Ref. 99 49 13  
B. 1 x 25 ml. Substrate. Ref. 99 67 43

### Working reagent

Reagents A and B are ready-to-use. If a Monoreagent procedure is preferred, then the reagents must be mixed in the ratio: 4 parts of A (Buffer solution) + 1 part of B (Liquid substrate).

2-Amino-2-methyl-1-propanol buffer pH 10.4 0.7 M  
p-Nitrophenylphosphate 12 mM  
HEDTA 1.55 mM  
Mg Acetate 1.5 mM  
Zinc Sulphate 0.95 mM  
Preservatives and stabilizers

### Storage and stability

The components of the kit, stored at 2-8°C, will remain stable until the expiration date stated on the label. The Monoreagent is stable for 3 weeks at 2-8°C and for 1 week at room temperature ( $\leq 25^\circ\text{C}$ ), when protected from the sunlight.

### Sample

Serum or heparinized plasma. Use samples free from hemolysis. Sera kept in the refrigerator at 2-8°C will remain stable for 7 days.

### Caution

The reagent contains Sodium azide at 0.09%. Handle with care. The disposal of the residues has to be made according to legal local regulations.

Monoreagent Procedure	ml
Working reagent	1.0
Sample	0.02

Mix, read the absorbance after 1 min. and start the stopwatch. Read again the absorbance after 1, 2, and 3 min.

Bireagent Procedure	ml
Buffer solution. (A)	1.0
Sample	0.1

Mix, incubate for approx. 1 minute and then add:

Substrate (B)	0.25
---------------	------

Mix, read the absorbance after 1 min. and start the stopwatch. Read again the absorbance after 1, 2, and 3 min.

### Reading

Wavelength: Hg 405 nm.  
Blank: Water.  
Cuvette: Thermostatzed 1 cm light-path.

### Calculations

Determine the  $\Delta E/\text{min.}$  For every reading and find the mean value.  $U/L = (\Delta E/\text{min.}) \times \text{Factor.}$

**Monoreagent**  
405 nm Factor: **2760**

**Bireagent**  
405 nm Factor: **3432**

### Normal values

Adults  
Women: 42 – 141 U/L  
Men: 53 – 128 U/L  
Children: Levels are very changeable. Up to 18 years old, its possible to find values between 42 – 383 U/L.  
In pregnancy women the Alkaline phosphatase value can be the double than not pregnancy women.

### Performance Characteristics

Linearity: up to  $\Delta E/\text{min.}$  values of 0.250. For higher values it is recommended to dilute the sample 1/10 with normal saline (NaCl 0.9%) and assay it once again.

Multiply the final result by 10. The analytical performance characteristics of the product depend both of the reagent and the reading system used, manual or automatic. The following data have been obtained manually.

Intraseres Variation Coefficient: 2.09 %.  
Interseres Variation Coefficient: 2.14 %.  
Recovery: 98.5 %.

Hemolysis and lipemia will interfere the assay. Anticoagulants such as EDTA, oxalate or citrate which chelate divalent cations should not be used since they would result in enzyme inhibition.

### Quality control

Seriscann Normal (Normal Control Serum) (Ref. 99 41 48) and Seriscann Anormal (Abnormal Control Serum) (Ref. 99 46 85).

### Autoanalyzers

Technical bulletins for different analyzers are available upon request.

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