FOSFORO INORGANICO U.V.

METODO DE FISKE – SUBBAROW

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INORGANIC PHOSPHORUS U.V.

FISKE – SUBBAROW METHOD

For "in vitro" determination of Phosphorus in serum, plasma or urine

Principle

The phosphate ion reacts with molybdate to produce phosphomolybdate which is finally reduced to a molybdenum blue, which is photometrically measured in the UV range.

Reagents

Kit 3 x 100 ml. (Ref. 99 76 09). Contents:

- A. 3 x 100 ml Reagent. Ref. 99 10 67
B. 1 x 5 ml Standard. Ref. 99 00 20

Aqueous solution of Phosphorus equivalent to 4 mg/dl (1.29 mmol/l).
Ready-to-use.

The reagent composition is as follows:

Ammonium molybdate	12 mM
Sulfuric acid	2.2 N
Preservatives and stabilizers	

Storage and stability

When stored at 2 -8° C, the reagent will remain stable until the expiration date stated on the label. Likewise, the standard will remain stable until the expiration date stated on the label, if kept at room temperature ($\leq 25^{\circ}\text{C}$).

Sample

Serum, heparinised plasma or urine.

Do not use plasma obtained with oxalate, citrate or EDTA, because can interfere the test.

When assaying urine samples, they shall be diluted 1/20 in deionized water prior to testing. Final result will be, therefore, multiplied by 20.

Serum phosphate remains stable one week when kept at 2-8° C. In urine and at an acid pH, phosphate will become stable for up to 3 months if the sample is kept at 2-8° C.

Caution

The reagent is corrosive. Thus, security pipettes shall be used.

On the other hand, the standard contains sodium azide (0.09 %) as a preservative.

The disposal of the residues has to be made according to legal local regulations.

Procedure	BL ml	ST ml	PR ml
Sample	--	--	0.01
Standard	--	0.01	--
Reagent	1.00	1.00	1.00

Mix well and let stand for 10 min.at room temperature (20-25°C).

Reading
Wavelength: 340 nm.
Blank: BL contents.
Colour stability: a minimum of 1 hour.

Calculations
 $\frac{\text{SA O.D.}}{\text{ST O.D.}} \times 4 = \text{mg of phosphorus/dl}$

SI Units
(mg/dl) x 0.3229 = mmol/L

Normal values
Serum:
Children: 4.0 - 7.0 mg/dl
Adults: 2.5 - 5.0 mg/dl
Urine: 0.3 - 1.0 g/24 hours.

Performance. Characteristics.

The reaction is linear up to 15 mg/dl (4.84 mmol/L). If higher concentrations are to be run, dilute the sample with deionized water (1+1). Multiply the final result by 2.

The analytical performance characteristics of the product depend both of the reagent and the reading system used, manual or automatic. The following data have been obtained manually.

Intraseres Variation Coefficient: 2.33%
Interseres Variation Coefficient: 2.78%
Recovery: 96.4 %.

Lipemic and hemolyzed sera will interfere with the assay.

Detergents currently used in the laboratory contain phosphates and other substances that may interfere with the reaction. Therefore it is recommended to rinse the glassware with diluted nitric acid and finally with deionized water.

It is recommended to use disposable plasticware to run the tests to avoid any possible undesirable contamination.

Autoanalyzers

Technical adaptations for different analyzers are available upon request.

Quality control

Seriscann Normal (Ref. 99 41 48) and Seriscann Anormal (Ref. 99 46 85).

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