gamma - GT LÍQUIDA

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
gamma - GT LIQUID

SZASZ MODIFIED METHOD

For "in vitro" determination of gamma-GT in serum or plasma

Principle

Reagents

Kit 1 x 50 ml. (Ref. 99 35 61). Contents:

- A. 1 x 40 ml. Buffer solution. Ref. 99 92 48
B. 1 x 10 ml. Liquid substrate. Ref. 99 38 43

Kit 1 x 125 ml. (Ref. 99 10 56). Contents:

- A. 1 x 100 ml. Buffer solution. Ref. 99 10 68
B. 1 x 25 ml. Liquid substrate. Ref. 99 03 93

Kit 1 x 250 ml. (Ref. 99 10 54). Contents:

- A. 2 x 100 ml. Buffer solution. Ref. 99 10 68
B. 1 x 50 ml. Liquid substrate. Ref. 99 10 82

Working reagent

Reagents A and B are ready-to-use. If a Monoreagent procedure is preferred, then the reagents must be mixed in the ratio: 4 parts of A (Buffer solution) + 1 part of B (Liquid substrate).

The concentrations in the reagent solution are:

Tris buffer pH 8.4	90 mM
L-gamma-glutamyl-3-carboxy-p-nitroanilide	6 mM
Glycyl-glycine	100 mM
Stabilizers and preservatives	

Storage and stability

The components of the kit, stored at 2-8°C, will remain stable until the expiration date stated on the label.
The Monoreagent is stable for 4 weeks at 2-8°C and for week at room temperature ($\leq 25^\circ \text{C}$), when protected from the sunlight.

Sample

Serum or plasma. Samples free from hemolysis should be used. The gamma-GT-asic activity will remain stable for 10 days if the sample is stored in the refrigerator at 2-8°C.

Caution

The substrate as well as the buffer solution contain sodium azide (0.09 %) as a preservative. The disposal of the residues has to be made according to legal local regulations.

Monoreagent procedure	ml
Working reagent	1.0
Incubate 2-3 min. at the desired assay temperature (25, 30 or 37°C).	
Sample	0.1
Mix, read the absorbance after 1 min. and start the stopwatch. Read again the absorbance after 1, 2 and 3 min.	
Bireagent procedure	ml
Enzymes sol. (A)	1.0
Sample	0.1
Mix, incubate for approx. 1 minute and then add:	
Substrate (B)	0.25
Mix, read the absorbance after 1 min. and start the stopwatch. Read again the absorbance after 1, 2 and 3 min.	

Reading

Wavelength: 405 nm.
Blank: Water.
Cuvette: Thermostated 1 cm light-path.

Calculations

Determine the $\Delta E/\text{min.}$ for every reading and find the mean value:

$$\text{U/L} = (\Delta E_{405 \text{ nm/min.}}) \times 1158 \text{ (Monoreagent)}$$
$$\text{U/L} = (\Delta E_{405 \text{ nm/min.}}) \times 1420 \text{ (Bireagent)}$$

Normal Values

Assay temp.	Men	Women
25 °C	6-28 U/L	4-18 U/L
30 °C	8-38 U/L	5-25 U/L
37 °C	11-50 U/L	7-32 U/L

Performance Characteristics

Linearity: Up to $\Delta E/\text{min.} = 0.200$. For higher values, dilute the sample 1/10 in saline (NaCl 0.9%) and assay once again. Multiply the final result by 10.

The analytical performance characteristics of the product depend both of the reagent and the reading system used, manual or automatic. The following data have been obtained manually.

Intraseres Variation Coefficient: 1.77%
Interseries Variation Coefficient: 2.29%
Recovery: 98.5 %.

Hemolyzed samples will give false results.
Fluoride, EDTA, citrate and oxalate inhibit the enzyme activity.

Quality control

Seriscann Normal (Normal Control Serum) (Ref.99 41 48) and Seriscann Anormal (Abnormal Control Serum) (Ref.99 46 85).

Autoanalyzers

Technical bulletins for different analyzers are available upon request.

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