GLUCOSA LIQUIDA

METODO GOD – POD

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
GLUCOSE LIQUID

GOD – POD METHOD

For “in vitro” determination of Glucose in serum or plasma

Principle

Reagents

Kit 1 x 100 ml (Ref. 99 82 25) Contents:

A. 1 x100 ml. Reagent. Ref. 99 82 84
B. 1 x 5 ml. Standard. Ref. 99 02 93

Kit 3 x 100 ml (Ref. 99 82 82) Contents:

A. 3 x100 ml. Reagent. Ref. 99 82 84
B. 1 x 5 ml. Standard. Ref. 99 02 93

Kit 4 x 250 ml (Ref. 99 86 60) Contents:

A. 4 x250 ml. Reagent. Ref. 99 01 68
B. 1 x 5 ml. Standard. Ref. 99 02 93

Working reagent

The reagent is ready-to-use.

Concentrations in the reagent solution are:

Phosphate buffer pH 6.8	100 mM
p-Hydroxybenzoic acid	39.5 mM
4-Aminoantipyrine	0.8 mM
Phenol	4.5 mM
Glucose Oxidase	≥ 18 kU/l
Peroxidase	≥ 1.1 kU/l
Preservatives and stabilizers	

Standard: Aqueous solution equivalent to 100 mg/dl (5.55 mmol/L).

Storage and stability

The components of the kit, stored at 2-8°C, will remain stable until the expiration date stated on the label.

Sample

Serum, plasma or C.S.F. Serum or plasma glucose (but not whole blood glucose, due to glycolytic processes) can be stored during 2 - 3 days, when refrigerated at 2-8°C.

Caution

The reagent contains Sodium azide at 0.09%. Handle with care. The disposal of the residues has to be made according to legal local regulations. Do not pipette directly from the bottle of reagent, for avoid undesired contaminations.

Performance Characteristics

The assay is linear up to 500 mg of Glucose/dl. For higher concentrations dilute the sample 1/2 with saline (NaCl 0.9%). Multiply the final result by 2. The analytical performance characteristics of the product depend both of the reagent and the reading system used, manual or automatic. The following data have been obtained manually.

Intraseres Variation Coefficient: 0.79%
Interseries Variation Coefficient: 1.33%
Recovery: 98.9 %.

Interferences:

Haemoglobin, higher than 200 mg/dl
Bilirubin, higher than 20 mg/dl;
Uric acid, higher than 20 mg/dl
Creatinine, higher than 15 mg/dl.

Interferences caused by the anticoagulants of current use such as Heparin, EDTA or Oxalate have not been described.

Quality control

Seriscann Normal (Normal Control Serum) (Ref. 99 41 48) Seriscann Abnormal (Abnormal Control Serum) (Ref. 99 46 85).

Autoanalyzers

Technical bulletins for different analyzers, are available upon request.

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Valores normales

Suero, plasma: 75 - 115 mg/dl.
L.C.R.: 40 - 80 mg/dl.

Procedure	BL	SA	ST
	ml	ml	ml
Standard	---	---	0.01
Sample	---	0.01	---
Working reagent	1.00	1.00	1.00

Mix well and incubate 5 - 10 min. at 37°C or 20-25 min. at 15 – 25°C.

Reading

Wavelength: Hg 546 nm; 505 nm.
Blank: the contents of BL.
Colour stability: a minimum of 1 hour, when protected from direct sunlight.

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