LDH LIQUIDA

METODO SFBC

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
LDH LIQUID

SFBC METHOD

For "in vitro" determination of Lactate Dehydrogenase in serum or plasma

Principle

Reagents

Kit 1 x 50 ml. (Ref. 99 18 18). Contents:

- A. 1 x 40 ml. Buffer solution. Ref. 99 13 74
B. 1 x 10 ml. Buffered NADH. Ref. 99 32 19

Kit 1 x 125 ml. (Ref. 99 00 35). Contents:

- A. 1 x 100 ml. Buffer solution. Ref. 99 00 40
B. 1 x 25 ml. Buffered NADH. Ref. 99 00 45

Working reagent

Reagents **A** and **B** are ready-to-use. If a **Monoreagent** procedure is preferred, then the reagents must be mixed in the ratio: 4 parts of **A** (Buffer solution) + 1 part of **B** (Buffered NADH).

The concentrations in the reagent solution are:

Tris-HCl buffer pH 7.2	80 mM
Sodium pyruvate	1.6 mM
NaCl	200 mM
NADH	0.20 mM
Stabilizers and preservatives	

Storage and stability

The components of the kit, stored at 2-8°C, will remain stable until the expiration date stated on the label. The **Monoreagent** is stable for 4 weeks at 2-8°C and for 1 week at room temperature ($\leq 25^\circ\text{C}$), when protected from the sunlight.

Sample

Serum, heparinized plasma. Samples free from hemolysis should be used. The enzyme in serum is stable for 2 days at 2-8°C. Frozen samples are rapidly inactivated.

	30°C / 37°C
Monoreagent procedure	ml
Working Reagent	1.0
Sample	0.02
Mix, read the absorbance after 1 min. and start the stopwatch. Read again the absorbance after 1, 2 and 3 min.	
Técnica Birreactivo	ml
Buffer solution (A)	1.0
Sample	0.02
Mix, incubate for approx. 1 min., then add:	
NADH (B)	0.2
Mix, read the absorbance after 1 min. and start the stopwatch. Read again the absorbance after 1, 2 and 3 min.	
Reading	
Wavelength: Hg 334 nm; 340 nm; Hg 365 nm.	
Blank: Water.	
Cuvette: Thermostatzated, 1 cm light-path.	
Calculations	
Determine the $\Delta E/\text{min.}$ for every reading and find the mean value.	
$U/L = (\Delta E/\text{min.}) \times \text{Factor.}$	

Caution

The reagent contains Sodium azide at 0.09%. Handle with care. The disposal of the residues has to be made according to legal local regulations.

Performances Characteristics

Linearity: Up to $\Delta E/\text{min.}$ of 0.150 (at 340 and 334 nm) or of 0.080 (at 365 nm). For higher values it is recommended to dilute the sample 1/10 with normal saline and assay it once again. Multiply the final result by 10. The analytical characteristics of the product depend both on the reagent and on the reading system used, manual or automatic. The following data have been obtained manually.

Intraseres Variation Coefficient: 1.90%
Interseres Variation Coefficient: 2.14%
Recovery: 98.8 %.

It is recommended to use heparin as anticoagulant. Others, like citrate, oxalate, or fluoride may interfere with the test. Hemolized samples will give false results.

When assaying high activity samples, a very low initial absorbance can be found, which is mainly due to the rapid consumption of the NADH at the early stage of the reaction. In such a case, dilute the sample with saline (NaCl 0.9 %) and assay it once again.

Quality control

Seriscann Normal (Normal Control Serum) (Ref. 99 41 48) and Seriscann Anormal (Abnormal Control Serum) (Ref. 99 46 85).

Autoanalyzers

Technical bulletins for different analyzers, are available upon request.

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