# LIPASA

## METODO COLORIMETRICO

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# LIPASE

## COLORIMETRIC METHOD

For the "in vitro" determination of Lipase in serum or plasma

### Principle

At alkaline pH the lipase dye substrate 1,2-O-dilauryl-rac-glycero-3-glutaric acid-(6-methyl-resorufin) ester is cleaved by the catalytic action of pancreatic lipase to form 1,2-O-dilauryl-rac-glycerol and unstable compound glutaric acid-(6-methyl-resorufin)-ester. In that alkaline solution, this compound degrades to glutaric acid and methyl-resorufin. The color intensity of the red dye formed is directly proportional to the lipase activity. This activity can be measured at 578 nm.

### Reagents

Kit 1 x 80 ml (Ref. 99 11 15). Contains:  
A. 1 x 50 ml. Buffer solution. Ref. 99 11 20  
Ready-to-use.  
B. 1 x 30 ml. Substrate. Ref. 99 12 15  
Ready-to-use.  
C. 1 x 1ml. Freeze-dried Standard. Ref. 99 13 07  
Rehydrate the vial with 1 ml. of deionized water.  
The concentration is stated on the label of the vial.

### Working reagent

The concentrations in the reagents solutions are:  
A. Buffer solution 50 mM  
Bicin buffer pH 8  $\geq 1$  mg/l  
Co-lipase 1,8 mM  
Sodium deoxycholate 12 mM  
Calcium chloride

B. Substrate  
Tartrate buffer pH 4 12 mM  
1,2-O-dilauryl-rac-glycero-3-glutaric acid-(6-methylresorufin) - ester 0.27 mM  
Taurodeoxy-cholate 9.0 mM  
Surfactants and preservatives

### Storage and stability

When stored at 2° - 8° C, the components of the kit will remain stable until the expiration date stated on the label.  
The standard, once rehydrated, is stable for 15 days at 2-8°C

### Sample

Serum or heparinized plasma. Use samples free from hemolysis. Do not use plasma with EDTA.  
The Lipase activity in serum is stable for 4 days at 2°-8° C and for 24 h at room temperature (20 – 25°C).

### Caution

The reagents contain Sodium azide at 0.09%. Handle with care. The disposal of the residues has to be made according to legal local regulations.

Procedure	BL	PR	ST
	ml	ml	ml
Water	0,01	--	--
Sample	--	0,01	--
Standard	--	--	0,01
<b>Reagent A</b>	1,0	1,0	1,0
Mix well and let stand for 5 min / 37°C. Add Reagent B.			
<b>Reagent B</b>	0,6	0,6	0,6
Mix well. After 60 seconds read (E1). After 90 seconds more, read again (E2).			

### Reading

Wavelength: 578 nm.  
Blank: Contents of BL.

### Calculation

$$\Delta E = E2 - E1$$

$$\frac{\Delta E \text{ PR} - \Delta E \text{ BL}}{\Delta E \text{ ST} - \Delta E \text{ BL}} \times \text{STD Conc.} = \text{U/L.}$$

To express the results in  $\mu\text{kat/L}$ , multiply the U/L by the factor 0.0167.

### Normal values

Serum, plasma: < 60 U/L.  
Nevertheless each laboratory should investigate its own reference values.  
To convert the color units in turbidimetric units, multiply the final result by 3.2.

### Performance Characteristics

Linearity: Up to 300 U/L. For higher concentrations, dilute the sample 1/5 with saline (NaCl 0.9 %) and assay once again. Multiply the final result by 5.  
The analytical performance characteristics of the product depend both of the reagent and the reading system used, manual or automatic. The following data have been obtained manually.

Intraseres Variation Coefficient: 2.28%  
Interseries Variation Coefficient: 2.64%  
Recovery: 97.6 %.

Hemolysis interferes with the assay.  
Concentrations of Hemoglobin higher than 400 mg/dl and of Bilirubin higher than 60 mg/dl can interfere with the assay.  
To perform this test, the use of disposable labware is highly recommended.

### Quality control

Seriscann Normal (Normal control serum) (Ref.99 41 48) and Seriscann Anormal (Abnormal control serum) (Ref. 99 46 85).

### Autoanalyzers

Technical bulletins for different automatic analyzers are available upon request.

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