# MAGNESIO

## METODO CON AZUL DE XILIDILO

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# MAGNESIUM

## XYLIDYL BLUE METHOD

For "in vitro" determination of Magnesium in serum, plasma C.S.F. and urine

### Principle

In an alkaline medium, the magnesium ions of the sample will produce a coloured complex with Xylidyl blue. Colour intensity is directly proportional to the magnesium ions' concentration present in the sample. Glycoetherdiamine-N,N,N'- tetraacetic acid (EGTA) performs as a chelating agent.

### Reagents

**Kit 2 x 100 ml (Ref. 99 29 85).** Contents:

- A. 2 x 100 ml. Xylidyl blue reagent. Ref. 99 35 10  
B. 1 x 5 ml Standard. Ref. 99 06 40

Aqueous solution of Magnesium equivalent to 4 mg/dl (1.645 mmol/l).  
Ready-to-use.

The reagent composition is as follows:

Xylidyl blue	0.1 mM
Tris buffer pH 11	0.3 mM
EGTA	50 µM
K <sub>2</sub> CO <sub>3</sub>	75 mM
Surfactants	2%
Preservatives and stabilizers	

### Storage and stability

When stored at room temperature (15 - 25° C), the components of this kit will remain stable until the expiration date stated on the label.

### Sample

Non-hemolyzed serum, heparinized plasma, C.S.F., and urine. Serum magnesium is stable for up to 5 – 7 days at 2-8° C.

To run a urine sample, this should be previously taken to an acid pH value (3-4) by adding some drops of HCl. Then, dilute with deionized water (1+4). Multiply the final result by 5.

### Caution

The reagent and Standard contains Sodium azide at 0.09%. Handle with care. The disposal of the residues has to be made according to legal local regulations.

Procedure	BL	ST	PR
Sample	ml	ml	ml
Standard	--	--	0,01
Reagent	1,00	1,00	1,00

Mix well and let stand 10 min. at room temperature (20-25°C).

**Reading**  
Wavelength: Hg 546 nm; 520 nm.  
Blank: Contents of BL.  
Colour stability: 1 hour.

**Calculations**  
SA Abs.  
----- x 4 = mg of magnesium / dl  
ST Abs.

**SI Units**  
(mg/dl) x 0.4113 = mmol/L

**Normal values**  
Serum,plasma: 1.9 - 2.5 mg/dl  
C.S.F.: 2.5 - 3.5 mg/dl  
Urine: 1.0 - 10 mg/dl  
24 h urine: 50 - 150 mg/24 h.

### Performance. Characteristics

Linearity: The reaction is linear up to 5 mg/dl. If higher concentrations are to be run, dilute the sample with deionized water (1+1). Multiply the final result by 2.

The analytical performance characteristics of the product depend both of the reagent and the reading system used, manual or automatic. The following data have been obtained manually.

Intraseres Variation Coefficient: 2.52%  
Interseres Variation Coefficient: 2.79%  
Recovery: 104.3 %.

Hemolyzed serum samples should be discarded.

It is recommended to use disposable plasticware to run the tests to avoid any possible undesirable contamination.

Lipaemic sera as well as a high bilirubin concentration (up to 20 mg/dl) do not interfere with the assay.

It is recommended to use disposable plasticware to run the tests to avoid any possible undesirable contamination.

### Autoanalyzers

Technical bulletins for different analyzers available upon request.

### Quality Control

Seriscann Normal (Normal Control Serum) (Ref. 99 41 48) and Seriscann Anormal (Abnormal Control Serum) (Ref. 99 46 85).

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