COLESTEROL - HDL DIRECTO

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Sugiuchi, H., et.al.(1995). Clin.Chem., 41, 717-723.

HDL - CHOLESTEROL DIRECT

COLOURIMETRIC METHOD **NEW SAMPLE VOLUME**

For "in vitro" determination of HDL - Cholesterol in serum or plasma

Principle

HDL-Cholesterol DIRECT is used for the determination of this cholesterol fraction, with no previous treatments, that is to say, without precipitation nor centrifugation.

The method is based on the property of a detergent to liberate the HDL fraction by solubilization, which then reacts with the chromogen, Cholesterol Esterase and Oxidase to give rise to a measurable colour at 600 nm. The use of a polyanion stabilizes the lipoproteins (VLDL, LDL and chylomicrons) by adsorption, and therefore cannot react with the enzymatic complex.

Reagents

Kit 1 x 40 ml (Ref. 99 48 58). Contents:

A. 1 x 30 ml Reagent (A). Ref. 99 48 60
B. 1 x 10 ml Reagent (B). Ref. 99 48 62
C. 1 x 1 ml Freeze-dried HDL-Calibrator. Ref. 99 03 72

Kit 1 x 80 ml (Ref. 99 80 58). Contents:

A.1 x 60 ml.Reagent (A). Ref. 99 80 60
B.1 x 20 ml.Reagent (B). Ref. 99 80 62
C.1 x 1 ml. Freeze-dried HDL-Calibrator. Ref. 99 03 72

Kit 1 x 400 ml (Ref. 99 57 88). Contents:

A.3 x 100 ml.Reagent (A). Ref. 99 33 47
B.1 x 100 ml.Reagent (B). Ref. 99 49 40
C. 1 x 1 ml Freeze-dried HDL-Calibrator. Ref. 99 03 72

Working reagent

Reagents (A) and (B) are ready-to-use.

The calibrator should be rehydrated with 1 ml of deionised water and let stand for 20 min prior to be used. Its concentration is stated on the vial label.

Composition of reagents is:

Reagent (A)
N,N-bis(4-sulphobutyl)-3-methylaniline 1.2 mM
Polyanion/polymer 0.8 mM
Stabilizers and preservatives

Reagent (B)
Cholesterol esterase ≥ 550 KU/L
Cholesterol oxidase ≥ 300 KU/L
POD ≥ 1500 KU/L
4-Aminoantipyrine 0.95 mM
Stabilizers and preservatives

| Procedure | SA | CAL |
|--|-----|-----|
| | µl | µl |
| Sample | 4.0 | -- |
| Calibrator | -- | 4.0 |
| Reagent A | 300 | 300 |
| Mix and incubate 5 min/37 °C. Read (E ₁) for the sample and the calibrator. | | |
| Reagent B | 100 | 100 |
| Mix and incubate 5 min/37 °C.Read, once again, (E ₂). | | |
| Reading Wavelength: 600 nm (546 nm - 640 nm). Blank: read against air. Cuvette: thermostized. 1 cm light-path. | | |

Storage and stability

When kept at 2 - 8 °C, the reagent will remain stable until the expiration date stated on the label. Reagents (A) and (B), once opened are stable 2 months when stored between 2° - 8° C and protected from direct sunlight. Once rehydrated, the calibrator will remain stable for 10 days when stored between 2° - 8° C.

Calculations

$$\frac{E_2 - E_1 \text{ sample}}{E_2 - E_1 \text{ calibrator}} \times [\text{CAL}]^* = \text{mg/dl}$$

S.I. Units

mg/dl x 0.0259 = mmol/L.

(*). See concentration on the vial's label.

Clinical interpretation

Low HDL-Cholesterol. Increases the risk for heart disease.

Men < 40 mg/dL
Women < 50 mg/dL

High HDL-Cholesterol. Reduces the risk for heart disease.

Men > 60 mg/dL
Women > 60 mg/dL

According to the recommendations of the European Atherosclerosis Society:

| | Lipid disorder |
|---|----------------|
| Colesterol < 200 mg/dl | NO |
| Triglycerides < 200 mg/dl | |
| Colesterol 200 - 300 mg/dl if HDL-Cholesterol < 35mg/dl | YES |
| Colesterol > 300 mg/dl | YES |
| Triglycerides > 200 mg/dl | |

As whit any diagnostic test procedure, results should be interpreted considering all other test results and clinical status of the patient.

Sample

Serum or heparinized plasma. In case the test is not to be carried out the same day, it is recommended to freeze the sample at - 15°C. Stored in this way, it will remain stable for up to 1 week.

Caution

The disposal of the residues has to be made according to legal local regulations.

Performance Characteristics

The analytical performance characteristics of the product depend both of the reagent and the reading system used, manual or automatic. The following data have been obtained using a Targa BT Analyzer.

The assay is linear up to 200 mg of Cholesterol/dl. For higher concentrations dilute the sample 1/2 with saline (NaCl 0.9%). Multiply the final result by 2.

Intraseres Variation Coefficient: 1.4%

Interseres Variation Coefficient: 2.2%

Bilirubin concentrations up to 40 mg/dl, Ascorbic acid up to 10 mM and Haemoglobin up to 1000 mg/dl and Triglyceride up to 1000 mg/dL will not interfere with the assay.

Quality control

Seriscann Normal. (Normal control serum, Ref. 99 41 48)

Seriscann Anormal. (Abnormal control serum, Ref. 99 46 85).

Autoanalyzers

Technical bulletins for different analyzers are available upon request.

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