COLESTEROL - HDL

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HDL - CHOLESTEROL

DEXTRAN SULPHATE - Mg (II) METHOD

For "in vitro" determination of HDL - Cholesterol in serum.

Principle

LDL and VLDL (low and very low density lipoproteins) are precipitated from serum by the action of a polysaccharide, in the presence of divalent cations. Then, high density lipoproteins cholesterol (HDL-Cholesterol) present in the supernatant, is determined.

Reagents

1 x 4 ml Precipitant solution. Ref. 99 38 85
Dropper for a minimum of 100 tests.
Ready-to-use.

Concentrations in the reagent solution are:

Dextran sulphate 10 g/L
Magnesium acetate 1 M
Stabilizers

Storage and stability

When kept at 2 - 8°C, the reagent will remain stable until the expiration date stated on the label.

Sample

Serum. The separation of the HDL fraction has to be carried out as soon as possible; otherwise, it is recommended to freeze the sample at - 15°C. Stored in this way, it will remain stable for up to 1 week.

Caution

The disposal of the residues has to be made according to legal local regulations.

Procedure

1. Precipitation reaction:

Sample 0,3 ml
Precipitant solution 1 drop

Mix and let stand for 15 min. at room temperature (20 - 25°C). Centrifuge at 2000 x g (1500 - 2300 rpm) / 15 min. or 10,000 x g (8000 - 12000 rpm) / 2 min.
Determine, in the supernatant, the concentration of cholesterol.

2.- Determination of Cholesterol. (1)

	BL ml	SA ml	ST ml
Supernatant	--	0.01	--
Standard	--	--	0.01
Working reagent	1.00	1.00	1.00

Mix well and let stand for 5 min at 37 °C or 10 min at room temperature

Reading

Wavelength: Hg 546 nm; 505 nm
Blank: the contents of BL
Colour stability: 1 hour (when protected from direct light)

Calculations

$\frac{\text{SA. Abs}}{\text{ST. Abs}} \times 200 \times 1,13 = \text{mg HDL-Cholesterol} / \text{dL}$

SI Units: (mg/dL) x 0.0259 = mmol/L

Remark: When the precipitant reagent is added to the sample, the latter is diluted by a factor of 1.13. Therefore, the cholesterol value should be multiplied by such a factor
(1) With QCA Cholesterol reagent

Clinical interpretation

Low HDL-Cholesterol. Increases the risk for heart disease.

Men < 40 mg/dL
Women < 50 mg/dL

High HDL-Cholesterol. Reduces the risk for heart disease

Men > 60 mg/dL
Women > 60 mg/dL

According to the recommendations of the European Atherosclerosis Society

		Lipid disorder
Cholesterol	< 200 mg/dL	NO
Triglycerides	< 200 mg/dL	NO
Cholesterol	200 - 300 mg/dL	SI
HDL Cholesterol	< 35 mg/dL	SI
Cholesterol	> 300 mg/dL	SI
Triglycerides	> 200 mg/dL	SI

Reagents necessary but not provided.

Reagent for Cholesterol determination and standard
It is advisable to use QCA reagents:

Kit 1 x 100 ml 99 52 82
Kit 3 x 100 ml 99 52 80
Kit 2 x 250 ml 99 50 12

Standard is included in the kits

Performance Characteristics

The assay is linear up to 700 mg of Cholesterol/dl. For higher concentrations dilute the sample 1/2 with saline (NaCl 0.9%). Multiply the final result by 2.

The analytical performance characteristics of the product depend both of the reagent and the reading system used, manual or automatic. The following data have been obtained manually.

Intraseres Variation Coefficient: 1.87%

Interseries Variation Coefficient: 2.22%

Recovery: 97.4 %.

Neither aged nor hemolyzed samples shall be assayed. Bilirubin concentrations higher than 9 mg/dl will interfere the precipitation. Sera with triglycerides concentration higher than 350 mg/dl will be diluted 1/2 in saline (NaCl 0.9%) prior to precipitation.

Remark

When the precipitant reagent is added to the sample, the latter is diluted by a factor of 1.13. Therefore, the cholesterol value should be multiplied by such a factor.

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