COLESTEROL – LDL DIRECTO

METODO COLORIMETRICO.

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Freedman, D., et al., (2004), Clin. Chem., **50**, 1189-1200.

Methods of Enzymatic Analysis, 3rd Edition, **VIII**, Edited by H. U. Bergmenyer. (1985), 154 - 160.

Bachorick, P.S., et al. (1995) Clin. Chem., **41**, Nº 10.

QUIMICA CLINICA APLICADA S.A.

Empresa certificada ISO 9001 / ISO 13485

A 7 Km 1081 – P.O. Box 20 - E43870 AMPOSTA / SPAIN

Tel. ++ 34 (977) 70. 62. 30 Fax ++ 34 (977) 70. 30. 40

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PRO4_REG9_LDLDI_03

LDL - CHOLESTEROL DIRECT

COLORIMETRIC METHOD

For "in vitro" determination of LDL - Cholesterol in serum

Principle

The LDL-Cholesterol reagent allows to measure the LDL Cholesterol fraction levels in serum or plasma directly, without the need for any sample pretreatment or centrifugation.

The method consists of two specific steps. In the first one, Chylomicron, VLDL and HDL are selectively eliminated by means of specific oxidation reactions with Cholesterol esterase and Cholesterol oxidase, and also by a specific degradation reaction with catalase, forming Cholestenone and H₂O₂.

In the second step the remainder of LDL Cholesterol can be spectrophotometrically measured when a color is formed as it is transformed (quinone complex) by the presence of specific surfactants.

Reagents

Kit 1 x 40 ml (Ref. 99 18 70). Contents:

A. 1 x 30 ml Reagent (A). (Ref. 99 18 72)

B. 1 x 10 ml Reagent (B). (Ref. 99 18 74)

C. 1 x 2 ml Freeze-dried LDL-Calibrator. (Ref. 99 18 78)

Rehydrate the vial with **2 ml** of deionized water. The concentration is stated on the label of the vial.

Working reagent

Reagents (A) and (B) are ready-to-use.

The calibrator, should be rehydrated with **2 ml** of deionized water and let stand for 20 min prior to be used. Its concentration is stated on the vial label.

The concentrations in the reagents are:

Reagent (A)

Good's buffer pH 7,0	100 mM
N-(2-hydroxy-3-sulphopropyl)-3,5-dimethoxyaniline	0.45 mM
Catalase	≥ 225 U/L
Cholesterol Esterase	≥ 600 kU/L
Cholesterol Oxidase	≥ 380 kU/L
Non reactive stabilizers	

Reagent (B)

Good's buffer pH 7,0	100 mM
4-Aminoantipirine	1.2 mM
Peroxidase	≥ 1200 kU/L
Non reactive stabilizers	

Storage and stability

When kept at 2 - 8°C, the reagent will remain stable until the expiration date stated on the label.

Reagents (A) and (B), once opened are stable for 2 months when stored at 2° - 8° C and protected from direct sunlight.

Once rehydrated, the calibrator will remain stable for 10 days when stored at 2° - 8° C.

Sample

Serum. The serum samples can be kept at 2 - 8°C for 1 day. Do not freeze.

Caution

The Reagent (B) contains sodium azide (0.09%) as a preservative. Handle with care. The disposal of the residues has to be made according to legal local regulations.

Procedure	BL µl	SA µl	CAL µl
Deionized water	5,00	---	---
Sample	---	5,00	---
Calibrator	---	---	5,00
Reagent A	300	300	300
Mix well and incubate at 37° C/5 min. Then add:			
Reagent B	100	100	100
Mix well and incubate at 37°C/ 5 min. Read (against the reagent's blank) the absorbance of the sample (Esa) and of the calibrator (Ecal).			

Calculations

(Esa)

_____ x conc. of calibrator = LDL- C (mg/dl).

(Ecal)

Reading

Wavelength: 600 nm.

Blank: the contents of BL.

Cuvette: 1 cm. light-path.

S. I. Units

mg/dl x 0.0259 = mmol/L.

Normal values

Low risk < 150 mg/dl.

Risk > 190 mg/dl.

It is advisable that each laboratory establishes its own reference values.

Performance Characteristics

Linearity: The assay is linear up to 1000 mg/dl.

The analytical performance characteristics of the product depend both of the reagent and the reading system used, manual or automatic. The following data have been obtained manually.

Intraseries Variation Coefficient: 1.31%

Interseries Variation Coefficient: 1.77%

Recovery: 99.1 %.

Concentrations of Ascorbic acid higher than 50 mg/dl, Hemoglobin higher than 500 mg/dl and of Bilirubin higher than 30 mg/dl can interfere with the assay.

Discard aged or hemolyzed sera.

Quality control

Seriscann Normal (Normal Control Serum)

(Ref. 99 41 48) and Seriscann Anormal (Abnormal Control Serum) (Ref. 99 46 85).

Autoanalyzers

Technical bulletins for different analyzers, are available upon request.

References

Freedman, D., et al., (2004), Clin. Chem., **50**, 1189-1200.

Methods of Enzymatic Analysis, 3rd Edition, **VIII**, Edited by H. U. Bergmeyer. (1985), 154 - 160.

Bachorick, P.S., et al. (1995) Clin. Chem., **41**, Nº 10.

QUÍMICA CLÍNICA APLICADA S.A.

ISO 9001 / ISO 13485 Certified Company

A 7 Km 1081 – P.O. Box 20 - E43870 AMPOSTA / SPAIN

Tel. ++ 34 (977) 70. 62. 30 Fax ++ 34 (977) 70. 30. 40

Revision: October 2009

PRO4_REG9_LDLDI_03