# COLESTEROL - LDL

## METODO CON SULFATO DE POLIVINILO

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# LDL - CHOLESTEROL

## POLYVINYL SULPHATE METHOD

For "in vitro" determination of LDL - Cholesterol in serum

### Principio

LDL-Cholesterol can be determined as the difference between Total Cholesterol and the Cholesterol content of the supernatant after precipitation of the LDL fraction by Polyvinyl sulphate (PVS) in the presence of Polyethylene-glycol monomethyl ether.

### Reagents

1 x 10 ml Precipitant solution. Ref. 99 06 10  
Dropper for 100 tests.  
Ready-to-use.

The concentrations in the reagent solution are:

Polyvinyl sulphate 0.7 g/L  
EDTA Na<sub>2</sub> 5.0 mM  
Polyethyleneglycol monomethyl ether 170 g/L  
Stabilizers

### Procedure

#### 1. Precipitation reaction

Precipitant solution 0.1 ml (3 drops)  
Sample 0.2 ml

Mix well and let stand for 15 min. aprox. at room temperature (20-25°C).

Centrifuge at 2,000 x g/15 min. (1500 - 2300 rpm) or at 10,000 x g/2 min. (8000 - 12000 rpm).

Determine the Cholesterol concentration in the supernatant.

#### 2. Cholesterol assay (1)

	BL	SA	ST
	ml	ml	ml
Supernatant	--	0.01	--
Standard	--	--	0.01
Working reagent	1.00	1.00	1.00

Mix well and let stand for 5 min at 37°C, or 10 min at room temperature.

### Reading

Wavelength: Hg 546 nm; 505 nm

Blank: The contents of BL

Colour stability: A minimum of 1 hour protected from the direct light.

### Calculations

#### a) Cholesterol in supernatant

$\frac{\text{Abs PR}}{\text{Abs ST}} \times 200 \times 1,5 = \text{mg supernatant Cholesterol/dl}$

#### b) LDL- Cholesterol

LDL- Cholesterol mg/dl = Total Cholesterol in the sample mg/dl - Supernatant Cholesterol mg/dl.

SI Units: mg/dl x 0.0259 = mmol/L

**Remark: When the precipitant reagent is added to the sample, the latter is diluted by a factor of 1.5. Therefore, the cholesterol value should be multiplied by such a factor.**

(1) With QCA Cholesterol reagent

### Storage and stability

When stored at 2 - 8°C the reagent will remain stable until the expiration date stated on the label.

### Sample

Serum. The serum samples can be kept at 2 - 8°C for 1 day. Do not freeze.

### Caution

The disposal of the residues has to be made according to legal local regulations.

### Clinical interpretation

Low risk for heart disease < 150 mg/dl.  
Risk for heart disease > 190 mg/dl.

### Performance Characteristics

The analytical performance characteristics of the product depend both of the reagent and the reading system used, manual or automatic. The following data have been obtained manually.

Intraseres Variation Coefficient: 2.87%  
Interseries Variation Coefficient: 2.96%  
Recovery: 96.4 %.

Falsely low values will sometimes occur if the assay is applied to sera with Triglyceride contents above 400 mg/dl. In such a case the sample shall be diluted 1/2 with saline (NaCl 0.9 %) prior to precipitation. Multiply the final result by 2.

Neither aged nor hemolyzed sera will be assayed

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