PROTEINAS TOTALES

METODO BIURET

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TOTAL PROTEINS

BIURET METHOD

For "in vitro" determination of Total Proteins in serum or plasma

Principle

At alkaline pH value proteins form a stable complex with Cu (II) ions, which is photometrically measured.

Reagents

Kit 3 x 100 ml. (Ref. 99 71 80). Contents:

A. 3 x 100 ml Biuret reagent.

Ref. 99 96 02

B. 1 x 5 ml Standard.

Ref. 99 02 46

The reagent composition is as follows:

NaOH	0.47 M
Potassium iodide	23.3 mM
Copper (II) sulphate	6.5 mM
Sodium-Potassium Tartrate	22.1 mM
Preservatives and stabilizers	

Standard. Aqueous solution of Proteins equivalent to 5 g/dl (50 g/L). Ready to use.

Storage and stability

When kept at room temperature (25°C), the reagent will remain stable until the expiration date stated on the label.

The Standard, however, should be stored at $2^{\circ} - 8^{\circ}\text{C}$.

It is therefore recommended, when using the 3x100 ml kit configuration, to take the standard out of the kit box and store it accordingly.

Sample

Serum or plasma, free from haemolysis. The sample will remain stable for up to 5 days, when kept at $2 - 8^{\circ}\text{C}$.

When using other corporal fluids as sample (amniotic liquid, urine, exudates, etc) it is advisable to take into account the margin values to test and the sensitivity of the reagent.

Caution

The reagent contains Sodium azide at 0.09%. Handle with care. The disposal of the residues has to be made according to legal local regulations.

Procedure	BL ml	ST ml	SA ml
Sample	--	--	0.02
Standard	--	0.02	--
Reagent	1.00	1.00	1.00

Mix well and let stand for 10 min. at room temperature ($20 - 25^{\circ}\text{C}$).

Reading
Wavelength: 540 nm.
Blank: BL contents.
Colour stability: a minimum of 3 hours.

Calculations
$$\frac{\text{SA O.D.}}{\text{ST O.D.}} \times 5 = \text{g of proteins/dl}$$

SI Units
(g/dl) x 10 = g/L

Normal values
6.6 - 8.7 g/dl

Performance Characteristics

Sensitivity: Detection limit 0,10 g/dl.

Reaction is linear up to values of 12 g/dl.

Concentrations higher than this value will require a prior dilution of the sample with saline (NaCl 0.9%).

Multiply by the dilution factor used.

The analytical performance characteristics of the product depend both of the reagent and the reading system used, manual or automatic. The following data have been obtained manually.

Intraserie Variation Coefficient: 0.85%

Interserie Variation Coefficient: 1.13%

Recovery: 98.7 %.

Quality control

Seriscann Normal (Normal Control Serum) (Ref. 99 41 48) and Seriscann Anormal (Abnormal Control Serum) (Ref. 99 46 85).

Autoanalyzers

Technical bulletins for different analyzers are available upon request.

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