TRIGLICERIDOS LIQUIDOS

METODO GPO

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
TRIGLYCERIDES LIQUID

GPO METHOD

For "in vitro" determination of Triglycerides in serum or plasma

Principle

Reagents

Kit 1 x 100 ml (Ref. 99 23 30). Contents

A. 1 x 100 ml. Reagent. Ref. 99 23 25
B. 1 x 5 ml. Standard. Ref. 99 03 17

Kit 3 x 100 ml (Ref. 99 23 20). Contents:

A. 3 x 100 ml. Reagent. Ref. 99 23 25
B. 1 x 5 ml. Standard. Ref. 99 03 17

Kit 2 x 250 ml (Ref. 99 30 80). Contents:

A. 2x 250 ml. Reagent. Ref. 99 01 53
B. 1 x 5 ml. Standard. Ref. 99 03 17

Working reagent

The reagent is ready-to-use.

Concentrations in the working reagent are:

Pipes buffer pH 6,8	50 mM
4-Chlorophenol	4.2 mM
4-aminoantipyrine	0.35 mM
ATP	2 mM
Magnesium aspartate	40 mM
Glycerol kinase	≥ 800 U/l
Glycerol-3-phosphate oxidase	≥ 2000 U/l
Peroxidase	≥ 500 U/l
Lipasas	≥ 9000 U/l
Non reactive Stabilizers	

Standard: Solution of Glicerol in water equivalent to 200 mg/dl (2.26 mmol/l). Ready-to-use

Storage and stability

The components of the kit, stored at 2-8°C, will remain stable until the expiration date stated on the label. Rehydrated vials are stable 60 days at 2 - 8°C and 10 days at room temperature (≤ 25°C), when protected from sunlight.

Procedure	BL	SA	ST
	ml	ml	ml
Sample	--	0.01	--
Standard	--	--	0.01
Working reagent	1.00	1.00	1.00

Mix well and let stand for 5 min. at 37° C, or 10 min at room temperature.

Reading
Wavelength: Hg 546 nm; 505 nm.
Blank: the contents of BL.
Colour stability: 1 hour. (when protected from direct sunlight).

Calculations
 $\frac{\text{SA O.D.}}{\text{ST O.D.}} \times 200 = \text{mg of Triglycerides/dl}$

SI Units
(mg/100 dl) x 0.0113 = mmol/L

Clinical interpretation

According to the recommendations of the European Atherosclerosis Society.

Lipid disorder	
Cholesterol < 200 mg/dl Triglycerides < 200 mg/dl	NO
Cholesterol 200 - 300 mg/dl if HDL cholesterol < 35 mg/dl	YES
Cholesterol > 300 mg/dl Triglycerides > 200 mg/dl	YES

Sample

Serum, or plasma with EDTA or heparin. Triglycerides are stable 4 days at 2-8°C, and up to three months at -20° C.

Caution

The reagent contains phenol derivatives. Handle with care. The disposal of the residues has to be made according to legal local regulations.

Performances Characteristics

The assay is linear up to 1000 mg/dl. Samples with a higher concentration shall be diluted 1/10 with NaCl 0.9% and assayed once again. Multiply the final result by 10.

The analytical performance characteristics of the product depend both of the reagent and the reading system used, manual or automatic. The following data have been obtained manually.

Intraseres Variation Coefficient: 0.89%
Interseries Variation Coefficient: 1.52%
Recovery: 98.5 %.

The assay is not interfered neither by hemoglobin (up to 150 mg/dl) nor bilirubin (up to 20 mg/dl).

Quality control

Seriscann Normal (Normal Control Serum) (Ref. 99 41 48) and Seriscann Anormal (Abnormal Control Serum) (Ref. 99 46 85).

Autoanalyzers

Technical bulletins for different analyzers, available upon request.

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