UREA ENZIMATICA

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
UREA ENZYMATIC

MODIFIED BERTHELOT – SEARCY METHOD

For "in vitro" determination of Urea in serum, plasma or urine

Principle

Reagents

Kit 5 x 100 ml (Ref.99 36 48). Contents:

A. 5 x 100 ml. Urease / Salicylate. Ref. 99 21 04

B. 1 x 15 ml Alkaline hypochlorite. Ref. 99 14 75

C. 1 x 5 ml Standard. Ref. 99 02 41

Aqueous solution of Urea equivalent to 40 mg / dl. (6.6 mmol/L).

Ready-to-use.

Working reagent

A. Dissolve the contents of the Urease/Salicylate vial with the volume of deionised water stated on the label.

The concentrations in the working reagent are:

Phosphate buffer pH 6.8	20 mM
Sodium salicylate	61 mM
Sodium nitroprusiate	3.4 mM
EDTA-Na2	1.34 mM
Urease	≥ 23 U/ml
Stabilizers	

B. Dilute the contents of the Alkaline hypochlorite vial up to 500 ml. of deionised water. Concentrations of reagent solution are:

Alkaline hypochlorite	7.5 mM
NaOH	160 mM

Storage and stability

The components of the kit, stored at 2-8°C, will remain stable until the expiration date stated on the label. Once the Urease/Salicylate vial has been dissolved, will remain stable for 3 weeks, if stored in amber bottle at 2-8° C. The Alkaline hypochlorite solution will remain stable for 8 months, if stored in the same way.

Sample

Serum, plasma and urine. Urea will remain stable in serum for at least 1 day at room temperature (≤ 25°C), 5 days at 2-8°C and 6 months when frozen (-20° C). In urine, urea will remain stable, when kept at 2-8°C, for 5 days, provided that the pH value be lower than 4.

If a urine sample is to be assayed, it should be previously diluted 1/100 with deionised water.

Multiply the final result by 100.

Caution

Reactivo B: In case of contact with the skin, mucose or eyes, wash thoroughly with water and ask the physician.

The reagent A contains Sodium azide at 0.09%. Handle with care. The disposal of the residues has to be made according to legal local regulations.

Procedure	BL	SA	ST
	ml	ml	ml
Standard	--	--	0.01
Sample	--	0.01	--
Reagent A	1.00	1.00	1.00
Mix and incubate 3 min. at 37°C. or 5 min. at room temperature (≤ 25°C).			
Reagent B	1.00	1.00	1.00
Mix and incubate again 3 min. at 37°C. or 5 min. at room temperature (≤ 25°C).			
Reading			
Wavelength: Hg 578 nm; 600 nm.			
Blank: The contents of the BL tube.			
Colour stability: 4 hours.			
Calculations			
$\frac{\text{SA O.D.}}{\text{ST O.D.}} \times 40 = \text{mg urea / dl}$			
S.I. Units			
(mg/dl) x 0.1665 = mmol/L			
Normal values			
Serum, plasma: 15 - 45 mg/dl.			
Urine: 20 - 35 g / L.			
Results as BUN (Blood Urea Nitrogen)			
mg/dL Urea x 0.467 = mg/dL BUN			

Performance Characteristics

Linearity: Up to 400 mg/dl of Urea. For higher values, dilute the sample 1/2 in deionised water and assay once again. Multiply the final result by 2.

The analytical performance characteristics of the product depend both of the reagent and the reading system used, manual or automatic. The following data have been obtained manually.

Intraserie Variation Coefficient: 1.66%
Interserie Variation Coefficient: 2.05%
Recovery: 97.9 %.

Any glassware contamination by ammonium salts or ammonia should be avoided. Serum samples should be free from haemolysis and turbidity. Fluoride as well as ammonium heparinate inhibit the reaction.

Quality control

Seriscann Normal (Normal Control Serum) (Ref. 99 41 48) and Seriscann Anormal (Abnormal Control Serum) (Ref. 99 46 85).

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