URIC ACID MONO SL

❖ Références/Referenes/ Referencias/ Referências:	Composition du coffret/ Kit composition/ Composición del kit/ Conteúdo da embalagem :
AUML-0420	R 6 x 50 mL
AUML-0500	R 6 x 100 mL
AUML-0250	R 12 x 20 mL
AUML-0427	R 6 x 50 mL + Std 1 x 5 mL
AUML-0497	R 1 x 100 mL + Std 1 x 5 mL
AUML-0507	R 6 x 100 mL + Std 1 x 5 mL
AUML-0707	R 4 x 250 mL + Std 1 x 5 mL

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English - EN

❖ INTENDED USE

ELITech Clinical Systems URIC ACID MONO SL is an *in vitro* diagnostic reagent intended for the quantitative determination of Uric Acid in human serum, plasma and urine samples.

CLINICAL SIGNIFICANCE ⁽¹⁻³⁾


Uric acid is the major product of the catabolism of endogenous and exogenous (dietary) purine nucleosides (adenosine and guanosine). This transformation mainly occurs in the liver. Approximately 75% of uric acid is eliminated by kidneys; the remainder is secreted into the gastrointestinal tract, where it is degraded by bacterial enzymes. Uric acid is not very soluble in water; urate crystals can occur in urines when the concentration is abnormally high. It can also happen in plasma, crystals then deposit essentially in joints, which induce intense inflammatory responses (gout). Some causes for increasing uric acid rate in serum are: increasing of purines synthesis, metabolic disorders (Lesch-Nyhan syndrome for example), nutritional troubles, increasing of nucleic acid turn-over in the case of proliferation of tumor cells, leukaemia, psoriasis, cytotoxic drugs, renal failures... Decreasing of uric acid rate in serum is more uncommon. It can occur in different cases: failure in renal elimination of uric acid (Fanconi syndrome), Hodgkin's disease for example. The quantitation of urinary uric acid is used to define the cause of hyperuricemia (excess purines or renal retention) and define appropriate treatment.

METHOD ⁽⁴⁾

Enzymatic - colorimetric.
 Trinder. End point.

PRINCIPLE ⁽⁴⁾

Enzymatic determination of uric acid according to the following reactions:

EHSPT = N-Ethyl-N-(2-Hydroxy-3-Sulfopropyl) *m*-Toluidine
 4-AAP = Amino-4-antipyrine

❖ COMPOSITION

Reagent: R		
Phosphate buffer, pH 7.0	100	mmol/L
EHSPT	0.72	mmol/L
Ferrocyanide	0.03	mmol/L
Amino-4-antipyrine	0.37	mmol/L
Uricase	≥ 150	U/L
Peroxidase	≥ 12 000	U/L
Sodium azide	< 0,1	%

Standard: Std (Ref : AUML-0427/0497/0507/0707)
 Uric acid 6 mg/dL

Sodium azide	<	357	µmol/L	%
		0.5		

❖ MATERIALS REQUIRED BUT NOT PROVIDED

- CALI-0550 ELICAL 2
- CONT-0060 ELITROL I
- CONT-0160 ELITROL II
- Normal saline solution (NaCl 9 g/L).
- General Laboratory equipment.
- Do not use materials that are not required as indicated above.

❖ WARNINGS AND PRECAUTIONS

- These *in vitro* diagnostic devices (reagent and standard) are for professional use only.
- The reagent R and standard contain sodium azide which may react with lead or copper plumbing to form potentially explosive metal azides. When disposing of these reagents always flush with copious amounts of water to prevent azide buildup.
- Take normal precautions and adhere to good laboratory practice.
- Use clean or single use laboratory equipment only to avoid contamination.
- The standard should be immediately and tightly capped to prevent contamination and evaporation.
- For more information, Safety Data Sheet (SDS) is available on request for professional user.

❖ STABILITIES

Store at 2-8 °C and protect from light. Do not freeze.
 Do not use after expiration dates indicated on the vial labels.

On board stability :
 The on-board stability is specific for each analyzer. (Refer to § PERFORMANCE DATA).

❖ PREPARATION

The reagent and standard are ready to use.

❖ PRODUCT DETERIORATION

- The product should be clear. Cloudiness would indicate deterioration.
- Do not use the product if there is visible evidence of biological, chemical or physical deterioration.
- Do not use the product if the damages of packaging might have an effect on the product performances (leakages, pierced vial).
- Do not use the product if there is visible evidence of biological, chemical or physical deterioration.
- Do not use the product if the damages of packaging might have an effect on the product performances (leakages, pierced vial).
- Do not use other specimens.

Warnings and precautions
 - According to Good Laboratory Practice, sampling should be performed prior to the administration of drugs. Sampling could lead to false results if performed during or immediately after the administration of some drugs.

- If unpreserved urine is received, add 0.1 mL of 12.5M NaOH to 10 mL of well-mixed urine; mix well. Warming at 60 °C to dissolve precipitate may be needed.

Storage and stability
 - Serum and heparinized plasma samples are stable 3 to 5 days if stored at 2-8 °C, 6 months at -20 °C.
 - Urines are stable 3 days at room temperature. Do not refrigerate urine samples.

REFERENCE VALUES ⁽⁵⁾

	Man	Woman	
Serum, plasma :	3.5 - 7.2	2.6 - 6.0	mg/dL
	208 - 428	155 - 357	µmol/L

Urine :	250 - 750	mg/24 h
	1.48 - 4.43	mmol/24 h
	16.7 - 50.0	mg/dL*
	0.99 - 2.97	mmol/L*

* for an urinary volume of 1.5 L per 24 hours

Note : The quoted range should serve as a guide only. It is recommended that each laboratory verifies this range or establishes a reference interval for the intended population.

❖ PROCEDURE

Manual Procedure	
Wavelength :	546 nm
Optical path :	1 cm
Sample/reagent ratio :	1:40
Temperature :	37 °C

Read against reagent blank.

	CALIBRATION	TEST
Reagent R	1 000 µL	1 000 µL
Standard/ Calibrator	25 µL	-
Sample	-	25 µL

Mix and read the absorbances (A) after an incubation of 5 minutes.

Automatic Procedure

These reagents may be used in several automatic analyzers. For ELITech Selectra Analyzers, validated applications are available on request. For Selectra TouchPro software, use the application included in the barcode available at the end of this insert.

❖ CALCULATION

A Sample x n n = Calibrator/ standard
 A Calibrator/ Standard concentration

Conversion factor : mg/dL x 59.48 = µmol/L
 mg/dL x 0.059 = mmol/L
 mg/dL x 10 = mg/L

Take the dilution factor into account for the calculation of uric acid concentration in urine.

❖ CALIBRATION

For the references AUML-0427/0497/0507/0707: For calibration, use either multiparametric calibrator ELICAL 2 or Uric Acid Standard 6 mg/dL. For the reference AUML-0250 0420/0500: For calibration, use multiparametric calibrator ELICAL 2.

Concentration values of Uric Acid Standard 6 mg/dL and multiparametric calibrator ELICAL 2 are traceable to the reference method ID-MS (Isotope Dilution - Mass Spectrometry).

Calibration frequency : The calibration is specific for each analyzer. (Refer to § PERFORMANCE DATA).

❖ QUALITY CONTROL

To check the accuracy of assays, control sera such as ELITROL I and ELITROL II should be used. These controls must be performed and validated before the patient samples are assayed. The control frequency must be at least once a day, after each calibration and should be adapted to Quality Control procedures of each laboratory and the regulatory requirements. Results should be within the defined ranges. If values fall outside of the defined ranges, each laboratory should take corrective measures. Quality control materials should be used in accordance with local guidelines.

WASTE MANAGEMENT

Disposal of all waste material should be in accordance with local, state and federal regulatory requirements.

PERFORMANCE DATA at 37 °C on ELITech Clinical Systems Selectra ProM Analyzers

❖- **Measuring range**
 Determined according to CLSI EP6-A protocol⁽⁶⁾.

a) **Serum/Plasma**
 The measuring range is from 1.50 to 25.00 mg/dL (89 to 1487 µmol/L).
 Samples having greater concentrations should be diluted 1:5 with NaCl 9 g/L solution and re-assayed. This procedure extends the measuring range up to 125.00 mg/dL (7436 µmol/L).

For users with Selectra TouchPro software, the «dilute» function performs the sample dilution automatically. Results take the dilution into account.

b) Urine

The measuring range is from 5.0 to 250.0 mg/dL (0.30 to 14.87 mmol/L).

Samples having greater concentrations should be diluted 1:5 with NaCl 9 g/L solution and re-assayed. This procedure extends the measuring range up to 800.0 mg/dL (47.59 mmol/L).

For users with Selectra TouchPro software, the «dilute» function performs the sample dilution automatically. Results take the dilution into account.

- **Limit of Detection (LoD) and Limit of Quantification (LoQ)**
 Determined according to CLSI EP17-A protocol⁽⁷⁾.

a) **Serum/Plasma**
 LoD = 0.09 mg/dL (5 µmol/L)
 LoQ = 1.00 mg/dL (59 µmol/L)

b) **Urine**
 LoD = 0.6 mg/dL (0.04 mmol/L)
 LoQ = 5.0 mg/dL (0.30 mmol/L)

- **Precision**
 Determined according to CLSI EP5-A2 protocol⁽⁸⁾.

a) **Serum/Plasma**

		Mean	Within-run	Total
	n	mg/dL	µmol/L	CV (%)
Low level	80	2.41	143	0.5 2.8
Medium level	80	4.95	294	0.7 2.3
High level	80	6.86	408	0.7 2.2

b) **Urine**

		Mean	Within-run	Total
	n	mg/dL	mmol/L	CV (%)
Low level	80	10.3	0.61	1.8 6.6
Medium level	80	23.9	1.42	1.1 3.8
High level	80	77.9	4.63	1.2 3.3

- Correlation

a) **Serum/Plasma**
 A comparative study has been performed between an ELITech Clinical Systems Selectra ProM Analyzer and another FDA-approved system equipment (Uricase method) on 100 human sera samples according to CLSI EP9-A2 protocol⁽⁹⁾. The values were between 1.55 and 23.94 mg/dL (92 and 1424 µmol/L). The parameters of the linear regressions are as follows:
 Correlation coefficient : (r) = 0.999
 Linear regression: y = 1.044 x - 0.04 mg/dL (2 µmol/L)

b) **Urine**
 A comparative study has been performed between an ELITech Clinical Systems Selectra ProM Analyzer and another FDA-approved system equipment (Uricase method) on 49 urine samples according to CLSI EP9-A2 protocol⁽¹¹⁾. The values were between 5.6 and 220.2 mg/dL (0.33 and 13.10 mmol/L). The parameters of the linear regressions are as follows:
 Correlation coefficient: (r) = 0.996
 Linear regression: y = 1.061 x + 0.1 mg/dL (0.01 mmol/L)

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- **Limitations/Interferences**
 - Do not use visibly turbid or hemolyzed samples.

- Do not report results outside of the usable range.

a) **Serum/ plasma**
 - Studies have been performed to determine the level of interference from different compounds according to CLSI EP7-A2 protocol⁽¹⁰⁾ and SFBC recommendations⁽¹¹⁾. Recovery is within ±10% of initial value at uric acid concentration of 2.52 mg/dL and 7.56 mg/dL.
Unconjugated bilirubin: No significant interference up to 30.0 mg/dL (513 µmol/L).
Conjugated bilirubin: No significant interference up to 14.8 mg/dL (252 µmol/L).
Hemoglobin: No significant interference up to 50 mg/dL.
Glucose: No significant interference up to 500 mg/dL (27.75 mmol/L).
Triglycerides: No significant interference up to 2095 mg/dL (23.67 mmol/L).
Turbidity: Interference occurs at all levels of Intralipid®
Ascorbic acid: Significant interference on samples containing ascorbic acid.
Methyl-dopa: No significant interference up to 1 mg/dL.
Calcium dobesilate: Induces falsely low results on individuals taking calcium dobesilate.

- In very rare cases, monoclonal gammopathies (multiple myeloma), in particular IgM type (Waldenström's macroglobulinemia) can cause unreliable results.⁽¹²⁾

- Results can be falsely lowered by significant levels in sample of NAC (N-Acetyl-Cysteine), NAPQI (metabolite of acetaminophene (paracetamol)) or metazolam.

- Many other substances and drugs may interfere. Some of them are listed in reviews published by Young.⁽¹³⁻¹⁴⁾

- The results of this assay should only be interpreted in conjunction with other diagnostic test results, clinical findings and the patient's medical history.

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**BIBLIOGRAPHIE/BIBLIOGRAPHY
BIBLIOGRAFÍA/BIBLIOGRAFIA**

1. First, M.R., *Renal Function. Clinical Chemistry: Theory, Analysis, Correlation*, 4th Ed., Kaplan, L.A., Pesce, A.J., Kazmierczak, S.C., (Mosby Inc. eds St Louis USA), (2003), 477 and appendix.
2. Newman, D. J., Price, C. P., *Non protein Nitrogen Metabolite*, Tietz Fundamentals of Clinical Chemistry, 5th Ed., Burtis, C.A. & Ashwood, E.R. (W.B. Saunders eds. Philadelphia USA), (2001), 414.
3. Tietz, N.W., *Clinical guide to laboratory tests*, 3rd Ed, (W.B. Saunders eds. Philadelphia USA), (1995), 624.
4. Fossati, P., et al., *Use of 3,5-dichloro-2-hydroxybenzenesul-fonic acid/4-aminophenazone chromogenic system in direct enzymatic assay of uric acid in serum and urine. Clin. Chem.*, (1980), **26**, 227.
5. Kaplan, L.A., Pesce, A.J., *Examination of Urine. Clinical Chemistry: Theory, Analysis, Correlation*, 4th Ed., Kaplan, L.A., Pesce, A.J., Kazmierczak, S.C., (Mosby Inc. eds St Louis USA), (2003), 1092.
6. *Evaluation of the Linearity of the Measurement of Quantitative Procedures: a Statistical Approach; Approved Guideline.* CLSI (NCCLS) document EP6-A (2003), **23** (16).
7. *Protocols for Determination of Limits of Detection and Limits of Quantification; Approved Guideline.* CLSI (NCCLS) document EP17-A (2004), **24** (34).
8. *Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline - Second Edition.* CLSI (NCCLS) document EP5-A2 (2004), **24** (25).
9. *Method Comparison and Bias estimation Using Patient Samples; Approved Guideline - Second Edition.* CLSI (NCCLS) document EP9-A2 (2002), **22** (19).
10. *Interference Testing in Clinical Chemistry; Approved Guideline - Second Edition.* CLSI (NCCLS) document EP7-A2 (2005), **25** (27).
11. Vassault A., et al., *Ann. Biol. Clin.*, (1986), 44, 686.
12. Berth, M. & Delanghe, J. *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature, Acta Clin Belg.*, (2004), **59**, 263.
13. Young, D. S., *Effects of preanalytical variables on clinical laboratory tests*, 2nd Ed., AACC Press, (1997).
14. Young, D. S., *Effects of drugs on clinical laboratory tests*, 4th Ed., AACC Press, (1995).

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- Symbols used are defined on ISO-15223-1 standard, except those presented below.

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- Only for ref. **AUML-0250**, used with Selectra TouchPro software.

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