FR

A partir du lot #18-xxxx et au-delà:

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English - EN

☛ INTENDED USE

For Total Bilirubin : ELITech Clinical Systems BILIRUBIN TOTAL 4+1 is an *in vitro* diagnostic reagent intended for the quantitative determination of total bilirubin in human serum and plasma samples on adults and children over 10 days.

For Direct Bilirubin : ELITech Clinical Systems BILIRUBIN DIRECT 4+1 is an *in vitro* diagnostic reagent intended for the quantitative determination of direct bilirubin in human serum and plasma samples.

CLINICAL SIGNIFICANCE (1-2)

Approximately 80-85 % of the bilirubin produced is derived from the heme moiety of the haemoglobin released from aging erythrocytes in the reticuloendothelial cells. Bilirubin, bound to albumin, is transported into the liver where it is rapidly conjugated with glucuronide to increase its solubility. Then it is excreted into biliary canaliculi, and hydrolyzed in the gastrointestinal tract.


Unconjugated bilirubin serum concentration increases in case of overproduction of bilirubin (acute or chronic hemolytic anaemia) and in case of disorders of bilirubin metabolism and transport defects (impaired uptake by liver cells: Gilbert's syndrome; defects in the conjugation reaction: Crigler-Najjar syndrome). Reduced excretion (hepatocellular damage: hepatitis, cirrhosis,...; Dubin-Johnson and Rotor syndrome) and obstruction to the flow of bile (most often produced by gallstones or by tumours) induce an important elevation of conjugated bilirubin and in a minor extent an increase of unconjugated bilirubin (conjugated hyperbilirubinemia).

METHOD (2)

Malloy-Evelyn modified. End point.

PRINCIPLE (2)

Sulphanilic acid reacts with sodium nitrite to form diazotized sulphanilic acid. In the presence of accelerator (cetrimide), conjugated and unconjugated bilirubin react with diazotized sulphanilic acid to form azobilirubin (Bilirubin Total 4+1). In the absence of accelerator, only conjugated bilirubin reacts (Bilirubin Direct 4+1). The increase of absorbance at 546 nm is proportional to bilirubin concentration.

COMPOSITION

BILIRUBIN TOTAL 4+1

Reagent 1: R1
 Sulphanilic acid 29 mmol/L
 Cetrimide 29 mmol/L

BILIRUBIN DIRECT 4+1

Reagent 1: R1
 Sulphanilic acid 29 mmol/L

BILIRUBIN TOTAL & DIRECT 4+1

Reagent 2: R2
 Sodium nitrite 11 mmol/L

☛ MATERIALS REQUIRED BUT NOT PROVIDED

- CALI-0550 ELICAL 2
 - CONT-0060 ELITROL I
 - CONT-0160 ELITROL II

- Normal saline solution (NaCl 9 g/L)
 - General Laboratory equipment.

- Biochemistry analyzer equipped with required filters (Refer to § PROCEDURE).

- Do not use materials that are not required as indicated above.

☛ WARNINGS AND PRECAUTIONS

- This *in vitro* diagnostic device is for professional use only.

- The reagents R1 contain sulphanilic acid . They may produce an allergic reaction.

- Take normal precautions and adhere to good laboratory practice.

- Use clean or single use laboratory equipment only to avoid contamination.

- Do not interchange reagent vials from different kits.

- For more information, Safety Data Sheet (SDS) is available on request for professional user.

STABILITIES

Store at 2-8 °C and protect from light. Do not freeze.

Do not use after expiration dates indicated on the vial labels.

On board stability:

The on-board stability is specific for each analyzer. (Refer to § PERFORMANCE DATA).

PREPARATION

The reagents are ready to use.

☛ PRODUCT DETERIORATION

- The reagent R1 of Bilirubin Total 4+1 can be slightly cloudy. It contains a detergent that can lead to the formation of foam in washing units of some equipments. These two characteristics are without consequences on the product performances.

- The reagent R1 of Bilirubin Direct 4+1 and the reagent R2 of Bilirubin Total & Direct 4+1 should be clear. Cloudiness would indicate deterioration.

- Do not use the product if there is visible evidence of biological, chemical or physical deterioration.

- Do not use the reagent if the damages of packaging might have an effect on the product performance (leakages, pierced vial).

SAMPLES (3)

Specimen

- Serum or lithium heparinized plasma.

- Do not use other specimens.

☛ Warnings and precautions

- **For Direct Bilirubin :** Do not use hemolysed samples.

- According to Good Laboratory Practice, sampling should be performed prior to the administration of drugs.

- Care must be taken to fill heparinized tubes according to the manufacturer's instructions. An insufficient filling may lead to erroneous results.

- Protect the samples from light before and during the analysis.

Storage and stability

If plasma and serum are protected from light, samples are stable 1 day (Total Bilirubin) or 2 days (Direct Bilirubin) at room temperature, 7 days at 2-8°C and 6 months at -20 °C.

REFERENCE VALUES (4)

Serum, plasma :

Total Bilirubin:

Adults and children over 10 days :
 0.2 - 1.2 mg/dL (3.4 - 21 µmol/L)

Direct Bilirubin:

< 0.2 mg/dL (3.4 µmol/L)

Note : *The quoted range should serve as a guide only. It is recommended that each laboratory verifies this range or establishes a reference interval for the intended population.*

☛ PROCEDURE

A) Total Bilirubin

Manual Procedure

Wavelength : 550 nm (λ. main) - 700 nm (λ. sub)

Optical path : 1 cm

Sample/reagent ratio : 1:20

Temperature: 37 °C

Read against distilled water.

	CALIBRATION	TEST
Reagent R1	800 µL	800 µL
Calibrator	50 µL	-
Sample	-	50 µL

Mix and read the absorbances (A) after an incubation of 5 minutes.

Calculate ΔA

$\Delta A = (A \lambda_{main}) - (A \lambda_{sub})$

	ΔA ₁	ΔA ₂
Then add :		
Reagent R2	200 µL	200 µL

Mix and read the absorbances (A) after an incubation of 5 minutes.

Calculate ΔA

$\Delta A = (A \lambda_{main}) - (A \lambda_{sub})$

	ΔA ₃	ΔA ₄
Automatic Procedure		

These reagents may be used in several automatic analyzers. For ELITech Selectra Analyzers, validated applications are available on request. For Selectra TouchPro software, use the application included in the barcode available at the end of this insert.

For Selectra ProXS users, a 700 nm additional filter is required.

B) Direct Bilirubin

Manual Procedure

Wavelength : 550 nm

Optical path : 1 cm

Sample/reagent ratio : 1:10

Temperature: 37 °C

Read against distilled water.

	CALIBRATION	TEST
Reagent R1	800 µL	800 µL
Calibrator	100 µL	-
Sample	-	100 µL

Mix and read the absorbances (A) after an incubation of 5 minutes.

	A ₁	A ₂
Then add :		
Reagent R2	200 µL	200 µL

Mix and read the absorbances (A) after an incubation of 50 seconds.

	A ₃	A ₄
Automatic Procedure		

Automatic Procedure

These reagents may be used in several automatic analyzers. For ELITech Selectra Analyzers, validated applications are available on request. For Selectra TouchPro software, use the application included in the barcode available at the end of this insert.

In the application, the offset must be set to : - 0.05 mg/dL (- 0.9 µmol/L).

☛ CALCULATION

A) Bilirubine Totale

$(\Delta A_4 - \Delta A_2) \times Fdil$
 x Calibrator concentration

$(\Delta A_3 - \Delta A_1) \times Fdil$

B) Direct Bilirubin

$(A_4 - A_2) \times Fdil$
 x Calibrator concentration

$(A_3 - A_1) \times Fdil$

Fdil = Dilution factor = (R1 volume + sample volume) / (R1 volume + R2 volume + sample volume)

Conversion factor: mg/dL x 17.1 = µmol/L

CALIBRATION

For calibration, multiparametric calibrator ELICAL 2 must be used. Its value is traceable to the reference material SRM 916a (of the National Institute of Standards and Technology).

Calibration frequency : The calibration is specific for each analyzer. (Refer to § PERFORMANCE DATA).

☛ QUALITY CONTROL

To check the accuracy of assays, control sera such as ELITROL I and ELITROL II should be used. These controls must be performed and validated before the patient samples are assayed. The control frequency must be at least once a day, after each calibration and should be adapted to Quality Control procedures of each laboratory and the regulatory requirements. Results should be within the defined ranges. If values fall outside of the defined ranges, each laboratory should take corrective measures. Quality control materials should be used in accordance with local guidelines.

WASTE MANAGEMENT

Disposal of all waste material should be in accordance with local, state and Federal regulatory requirements.

PERFORMANCE DATA at 37 °C on ELITech Clinical Systems Selectra ProM Analyzers

A) Total Bilirubin

☛ - Measuring range

Determined according to CLSI EP6-A protocol (5), the measuring range is from 0.25 mg/dL to 25.00 mg/dL (4.3 to 427.6 µmol/L).

Samples having greater concentrations should be diluted 1:5 with NaCl 9 g/L solution and re-assayed. This procedure extends the measuring range up to 60.00 mg/dL (1026.3 µmol/L).

For Selectra TouchPro users, the «dilute» function performs the sample dilution automatically. Results take the dilution into account.

- Limit of Detection (LoD) and Limit of Quantification (LoQ)

Determined according to CLSI EP17-A protocol (6),

LoD = 0.04 mg/dL (0.7 µmol/L)
 LoQ = 0.15 mg/dL (2.6 µmol/L)

- Precision

Determined according to CLSI EP5-A2 protocol (7).

		Mean		Within-run	Total
	n	mg/dL	µmol/L	CV (%)	
Level 1	80	1.15	19.7	1.8	5.0
Level 2	80	4.08	69.8	0.4	3.1
Level 3	80	14.61	249.9	0.5	2.9

- Correlation

A comparative study has been performed between an ELITech Clinical Systems Selectra ProM Analyzer and another FDA-Approved system equipment (DCA method) on 100 human serum samples according to CLSI EP9-A2 protocol (8). The sample concentrations were between 0.32 and 23.02 mg/dL (5.5 and 393.7 µmol/L). The parameters of the linear regressions are as follows:

Correlation coefficient: (r) = 0.999
 Linear regression: y = 0.948 x - 0.11 mg/dL (1.9 µmol/L)

- Limitations/Interferences

- Do not report results outside of the usable range.

- Studies have been performed to determine the level of interference from different compounds according to CLSI EP7-A2 protocol (9). Recovery within ± 15 % of initial value at total bilirubin concentration of 1.00 mg/dL and ± 10 % of initial value at total bilirubin concentration of 15.00 mg/dL.

Triglycérides: No significant interference up to 2100 mg/dL (23.73 mmol/L).

Hémoglobine: No significant interference up to 500 mg/dL.

Acetaminophen: No significant interference up to 30 mg/dL.

Ascorbic acid: No significant interference up to 4 mg/dL.

Acetylsalicylic acid: No significant interference up to 200 mg/dL.

- In very rare cases, monoclonal gammopathies (multiple myeloma), in particular IgM type (Waldenström's macroglobulinemia) can cause unreliable results. (10)

- Many other substances and drugs may interfere. Some of them are listed in reviews published by Young. (11-12)

- The results of this assay should only be interpreted in conjunction with other diagnostic test results, clinical findings and the patient's medical history.

- On board stability/Calibration frequency
On Board Stability: 28 days
Calibration frequency: 28 days

Recalibrate when reagent lots change, when quality control results fall outside the established range, and after a maintenance operation.

☛ These performances have been obtained using ELITech Selectra ProM analyzer. Results may vary if a different instrument or a manual procedure is used. The performances of applications not validated by ELITech are not warranted and must be defined by the user.

- The results of this assay should only be interpreted in conjunction with other diagnostic test results, clinical findings and the patient's medical history.

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BIBLIOGRAPHIE/BIBLIOGRAPHY

BIBLIOGRAFIA/BIBLIOGRAFIA

- Higgins, T., et al., Hemoglobin, Iron, and Bilirubin, Tietz Fundamentals of Clinical Chemistry, 6th Ed., Burtis, C.A. & Ashwood, E.R. (W.B. Saunders eds, Philadelphia USA), (2008), 509.
- Dufour, D.R. The liver: Function and chemical Pathology, Clinical Chemistry: Theory, Analysis, Correlation, 5th Ed., Kaplan, L.A., Pesce, A.J., (Mosby Inc.), (2010), 586 and appendix.
- Guder, W.G., et al., Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. (2002). WHO/DIL/LAB/99.1 Rev.2
- Wu, H.B., General Clinical Tests. Tietz Clinical guide to laboratory tests, 4th Ed., (W.B. Saunders eds, Philadelphia USA), (2006), 172.
- Evaluation of the Linearity of the Measurement of Quantitative Procedures: a Statistical Approach; Approved Guideline. CLSI(NCCLS) document EP6-A (2003), 23 (16).
- Protocols for Determination of Limits of Detection and Limits of Quantification; Approved Guideline. CLSI (NCCLS) document EP17-A (2004), 24 (34).
- Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline - Second Edition. CLSI (NCCLS) document EP5-A2 (2004), 24 (25) 8.
- Method Comparison and Bias estimation Using Patient Samples; Approved Guideline - Second Edition. CLSI (NCCLS) document EP9-A2 (2002), 22 (19).
- Interference Testing in Clinical Chemistry; Approved Guideline - Second Edition. CLSI (NCCLS) document EP7-A2 (2005), 25 (27).
- Berth, M. & Delanghe, J. Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature, Acta Clin Belg., (2004), 59, 263.
- Young, D.S., Effects of preanalytical variables on clinical laboratory tests, 2nd edition, AACCC Press (1997).
- Young D.S., Effects of drugs on clinical laboratory tests, 4th edition, AACCC Press (1995).

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

English - EN

INTENDED USE

For Total Bilirubin : ELITech Clinical Systems BILIRUBIN TOTAL 4+1 is intended for the quantitative *in vitro* diagnostic determination of total bilirubin in human serum and plasma on adults and children over 10 days.

For Direct Bilirubin : ELITech Clinical Systems BILIRUBIN DIRECT 4+1 is intended for the quantitative *in vitro* diagnostic determination of direct bilirubin in human serum and plasma.

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Unconjugated bilirubin serum concentration increases in case of overproduction of bilirubin (acute or chronic hemolytic anaemia) and in case of disorders of bilirubin metabolism and transport defects (impaired uptake by liver cells: Gilbert's syndrome; defects in the conjugation reaction: Crigler-Najjar syndrome).

Reduced excretion (hepatocellular damage: hepatitis, cirrhosis...; Dubin-Johnson and Rotor syndrome) and obstruction to the flow of bile (most often produced by gallstones or by tumours) induce an important elevation of conjugated bilirubin and in a minor extent an increase of unconjugated bilirubin (conjugated hyperbilirubinemia).

METHOD ⁽³⁾

Malloy-Evelyn modified. End point.

PRINCIPLE ⁽³⁾
 Sulfanilic acid reacts with sodium nitrite to form diazotized sulfanilic acid. In the presence of accelerator (cetrimide), conjugated and unconjugated bilirubin react with diazotized sulfanilic acid to form azobilirubin (Bilirubin Total 4+1). In the absence of accelerator, only conjugated bilirubin reacts (Bilirubin Direct 4+1). The increase of absorbance at 550 nm is proportional to bilirubin concentration.

Sulfanilic acid + NaNO₂ → Diazotized sulfanilic acid

Bilirubin + Diazotized sulfanilic acid → Azobilirubin

REAGENTS COMPOSITION

BILIRUBIN TOTAL 4+1

Reagent 1: R1
 Sulfanilic acid 29 mmol/L
 Hydrochloric acid 67 mmol/L
 Cetrimide 37 mmol/L

Reagent 2: R2
 Sodium nitrite 5,8 mmol/L

MATERIALS REQUIRED BUT NOT PROVIDED

- CALI-0550 ELICAL 2 4 x 3 mL
 - CONT-0060 ELITROL I 10 x 5 mL
 - CONT-0160 ELITROL II 10 x 5 mL
 - Normal saline solution (NaCl 9 g/L)
 - General Laboratory equipment.
 - Do not use materials that are not required as indicated above.

WARNINGS AND PRECAUTIONS

- These reagents are for professional *in vitro* diagnostic use only.

- Bilirubin Total reagent R1 is classified as hazardous :

WARNING : Causes skin irritation and serious eye irritation.

Wear protective gloves/protective clothing/eye protection/face protection. IF IN EYES : Rinse cautiously with water for several minutes.

Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists: Get medical advice/attention. IF ON SKIN: Wash with plenty of soap and water.

- The reagents R1 contain sulphanic acid . They may produce an allergic reaction.

- For more information, refer to the Safety Data Sheet (SDS).

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BILIRUBIN TOTAL & DIRECT 4+1

Reagent 2: R2
 Sodium nitrite 5,8 mmol/L

MATERIALS REQUIRED BUT NOT PROVIDED

- CALI-0550 ELICAL 2 4 x 3 mL
 - CONT-0060 ELITROL I 10 x 5 mL
 - CONT-0160 ELITROL II 10 x 5 mL
 - Normal saline solution (NaCl 9 g/L)
 - General Laboratory equipment.
 - Do not use materials that are not required as indicated above.

WARNINGS AND PRECAUTIONS

- These reagents are for professional *in vitro* diagnostic use only.

- Bilirubin Total reagent R1 is classified as hazardous :

WARNING : Causes skin irritation and serious eye irritation.

Wear protective gloves/protective clothing/eye protection/face protection. IF IN EYES : Rinse cautiously with water for several minutes.

Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists: Get medical advice/attention. IF ON SKIN: Wash with plenty of soap and water.

- The reagents R1 contain sulphanic acid . They may produce an allergic reaction.

- For more information, refer to the Safety Data Sheet (SDS).

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On board stability: The on-board stability is specific for each analyzer. (Refer to § PERFORMANCE DATA).

- The reagents R1 contain sulphanic acid . They may produce an allergic reaction.

- For more information, refer to the Safety Data Sheet (SDS).

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PROCEDURE

For ELITech Clinical Systems Selectra Analyzers: applications are available on request.

Wavelength : 546 nm

Temperature: 37 °C

Read against reagent blank.

Mix and read the absorbance (A1) after a 4 minutes and 40 seconds of incubation.

Mix and read the absorbance (A2) after a 6 minutes and 30 seconds of incubation.

With Selectra TouchPro software, use the application included in the barcode available at the end of this insert.

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Quality control : To ensure adequate quality, control sera such as ELITROL I (normal control) and ELITROL II (abnormal control) should be used. These controls must be performed and validated before the patient samples are assayed. The control frequency must be at least once a day, after each calibration and should be adapted to Quality Control procedures of each laboratory and the regulatory requirements. Results should be within the defined ranges. If values fall outside of the defined ranges, each laboratory should take corrective measures. Quality control materials should be used in accordance with local guidelines.

Waste management : Disposal of all waste material should be in accordance with local and legal requirements.

Performance data at 37 °C on ELITech Clinical Systems Selectra ProM Analyzers

A) Total Bilirubin

- **Measuring range**
 Determined according to CLSI⁽⁶⁾ EP6-A protocol, the measuring range is from 0.18 mg/dL to 6.82 mg/dL (3.1 to 116.7 µmol/L). Samples exceeding 6.82 mg/dL should be diluted 1:5 with NaCl 9 g/L solution (normal saline) and re-assayed. Use of this procedure extends the measuring range to 6.82 to 34.00 mg/dL (116.7 to 581.5 µmol/L).

- **Limit of Detection (LoD) and Limit of Quantification (LoQ)**
 Determined according to CLSI⁽⁷⁾ EP17-A protocol, LoD = 0.06 mg/dL (1.0 µmol/L)
 LoQ = 0.17 mg/dL (2.9 µmol/L)

- **Precision**
 Determined according to CLSI⁽⁸⁾ EP5-A2 protocol⁽⁸⁾.

B) Direct Bilirubin

- **Measuring range**
 Determined according to CLSI⁽⁶⁾ EP6-A protocol, the measuring range is from 0.18 mg/dL to 6.82 mg/dL (3.1 to 116.7 µmol/L). Samples exceeding 6.82 mg/dL should be diluted 1:5 with NaCl 9 g/L solution (normal saline) and re-assayed. Use of this procedure extends the measuring range to 6.82 to 34.00 mg/dL (116.7 to 581.5 µmol/L).

- **Limit of Detection (LoD) and Limit of Quantification (LoQ)**
 Determined according to CLSI⁽⁷⁾ EP17-A protocol, LoD = 0.06 mg/dL (1.0 µmol/L)
 LoQ = 0.17 mg/dL (2.9 µmol/L)

- **Precision**
 Determined according to CLSI⁽⁸⁾ EP5-A2 protocol⁽⁸⁾.

Correlation

A comparative study has been performed between an ELITech Clinical Systems Selectra ProM Analyzer and another FDA-Approved system equipment (DCA method) on 100 human serum samples according to CLSI EP9-A2 protocol⁽⁹⁾.

The sample concentrations were between 0.25 and 22.00 mg/dL (4.3 and 376.3 µmol/L).

The parameters of the linear regressions are as follows:
 Correlation coefficient: (r) = 0.999
 Linear regression: y = 0.924 x + 0.02 mg/dL (0.3 µmol/L)

Limitations/Interferences

- Ascorbic acid concentration above 2 mg/dL cause falsely elevated total bilirubin results.

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- Studies have been performed to determine the level of interference from different compounds according to CLSI EP7-A2 protocol⁽¹⁰⁾. Recovery within ± 10 % of initial value at total bilirubin concentration of 1.00 mg/dL and 15.00 mg/dL.

Triglycérides: No significant interference up to 2779 mg/dL (31.40 mmol/L).

Hemoglobin: No significant interference up to 500 mg/dL.

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- In very rare cases, monoclonal gammopathies (multiple myeloma), in particular IgM type (Waldenström's macroglobulinemia) can cause unreliable results.⁽¹¹⁾

- Many other substances and drugs may interfere. Some of them are listed in reviews published by Young.⁽¹²⁻¹³⁾

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- On board stability/Calibration frequency
On Board Stability: 28 days
Calibration frequency: 28 days
 Recalibrate when reagent lots change, when quality control results fall outside the established range, and after a maintenance operation.

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BIBLIOGRAPHIE/BIBLIOGRAPHY
BIBLIOGRAFIA/BIBLIOGRAFIA

- Nuttall, K.L., Klee, G.G., *Analytes of haemoglobin metabolism-Porphyrins, Iron and Bilirubin*, Tietz *Fundamentals of Clinical Chemistry*, 5th Ed., Burtis, C.A. & Ashwood, E.R. (W.B. Saunders eds. Philadelphia USA), (2001), 584.
- Schreiber, W.E., *Iron, Porphyrin and Bilirubin metabolism*, *Clinical Chemistry: Theory, Analysis, Correlation*, 4th Ed., Kaplan, L.A., Pesce, A.J., Kazmierczak, S.C., (Mosby Inc. eds St Louis USA), (2003), 657.
- Sherwin, J.E., Thompson, C., *Liver function*. *Clinical Chemistry: Theory, Analysis, Correlation*, 4th Ed., Kaplan L.A., Pesce, A.J., Kazmierczak, S.C., (Mosby Inc. eds St Louis USA), (2003), 493 and appendix.
- Guder, W.G., et al., *Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples*. (2002). WHO/DIL/LAB/99.1 Rev.2
- Wu, H.B., *General Clinical Tests. Tietz Clinical guide to laboratory tests*, 4th Ed., (W.B. Saunders eds. Philadelphia USA), (2006), 634-639.
- Evaluation of the Linearity of the Measurement of Quantitative Procedures: a Statistical Approach; Approved Guideline*. CLSI(NCCLS) document EP6-A (2003), 23 (16).
- Protocols for Determination of Limits of Detection and Limits of Quantification; Approved Guideline*. CLSI (NCCLS) document EP17-A (2004), 24 (34).
- Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline - Second Edition*. CLSI (NCCLS) document EP5-A2 (2004), 24 (25).8.
- Method Comparison and Bias estimation Using Patient Samples; Approved Guideline - Second Edition*. CLSI (NCCLS) document EP9-A2 (2002), 22 (19).
- Interference Testing in Clinical Chemistry ; Approved Guideline - Second Edition*. CLSI (NCCLS) document EP7-A2 (2005), 25 (27).
- Berth, M. & Delanghe, J. *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature*, *Acta Clin Belg.*, (2004), 59, 263.
- Young, D.S., *Effects of preanalytical variables on clinical laboratory tests*. 2nd edition, AACC Press (1997).
- Young D.S., *Effects of drugs on clinical laboratory tests*. 4th edition, AACC Press (1995).

SYMBOLES/SYMBOLS/ SIMBOLOS/SIMBOLOS

	Dispositif médical de diagnostic <i>in vitro</i> In vitro diagnostic medical device Dispositivo médico para diagnóstico <i>in vitro</i> Dispositivo médico de diagnóstico <i>in vitro</i>
	Consulter la notice d'utilisation Consult instruction for use Consulte el manual Consultar o manual de instruções
	Fabricant Manufacturer Fabricante Fabricante
	Limites de température Temperature limitation Límites de temperatura Limites de temperatura
	Numéro de lot Batch code Número de lote Número de lote
	Date d'expiration Use by Fecha de caducidad Prazo de validade
	Numéro de catalogue Catalogue number Número de catálogo Número de catálogo
	Contient Content Contiene Conteúdo
	Réactif 1 Reagent 1 Reactivo 1 Reagente 1
	Réactif 2 Reagent 2 Reactivo 2 Reagente 2
	Conformité Européenne European Conformity Conformidad Europea Conformidade Europeia

Bilirubin Total 0
220 PIT-BITD

Bilirubin Direct 0
200 PIT-BITD