**BIBLIOGRAPHIE/BIBLIOGRAPHY
BIBLIOGRAFIA/BIBLIOGRAFIA**

- Henderson, A.R., Donald W.M., *Enzymes, Tietz Fundamentals of Clinical Chemistry*, 5th Ed., Burtis, C.A. & Ashwood, E.R. (W.B. Saunders eds. Philadelphia USA), (2001), 352.
- Sanhai, W.R., Christenson, R.H., *Cardiac and muscle disease. Clinical Chemistry: Theory, Analysis, Correlation*, 4th Ed., Kaplan, L.A., Pesce, A.J., Kazmierczak, S.C., (Mosby Inc. eds St Louis USA), (2003), 566 and appendix.
- Neumeier, D., *et al., Clin Chim Acta.*, (1976), **73**, 445
- Schumann, G., *et al., Clin Chem Lab Med.*, (2002), **40**, 635.
- Klein, G., *et al., Clin Chem.*, (2001), **47**, Suppl. A30.
- Tietz N.W., *Clinical guide to laboratory tests*, 3rd Ed., W.B. Saunders eds. Philadelphia USA, (1995), 180.
- Vassault, A., *et al., Ann Biol Clin.*, (1986), **44**, 686.
- Vassault, A., *et al., Ann Biol Clin.*, (1999), **57**, 685.
- Berth, M. & Delanghe, J. *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature, Acta Clin Belg.*, (2004), **59**, 263.
- Young, D.S., *Effects of preanalytical variables on clinical laboratory tests*, 2nd Ed., AACCPress, (1997).
- Young, D.S., *Effects of drugs on clinical laboratory tests*, 4th Ed., AACCPress, (1995).

☛ **SYMBOLS/SYMBOLS/ SIMBOLOS/SIMBOLOS**

- Les symboles utilisés sont décrits dans la norme ISO-15223-1 hormis ceux présentés ci-dessous.
- Symbols used are defined on ISO-15223-1 standard, except those presented below.
- Los símbolos utilizados son descritos en la norma ISO-15223-1 a la excepción de los presentados a continuación.
- Os símbolos utilizados são definidos na norma ISO-15223-1, exceto os apresentados abaixo.

CONT	Contient Content Contiene Conteúdo
R1	Réactif R1 Reagent R1 Reactivo R1 Reagente R1
R2	Réactif R2 Reagent R2 Reactivo R2 Reagente R2
CE	Conformité Européenne European Conformity Conformidad Europea Conformidade Europeia

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English - EN

☛ **INTENDED USE**

ELITech Clinical Systems CK-MB SL is an *in vitro* diagnostic reagent intended for the quantitative determination of CK-MB form of creatine kinase (CK) in human serum samples.

CLINICAL SIGNIFICANCE ⁽¹⁻²⁾

Creatine kinase exists in 3 cytoplasmic forms: CK-MB (in cardiac muscle only), CK-MM (in striated and cardiac muscle) and CK BB (especially in brain). Determination of CK in serum is used for the diagnosis and the follow-up of cardiac muscle damages. In myocardial infarction, the rates of total CK and CK-MB increase quickly until reaching a peak 10-24h after the onset of the infarction. The levels return to normal within 3 to 4 days. CK-MB levels higher than normal can also be observed after muscle damages.

METHOD ⁽³⁻⁴⁾

Immuno-inhibition, IFCC Method. Kinetic. UV.

PRINCIPLE ⁽³⁻⁴⁾

CK-MB SL reagent contains an antibody inhibiting specifically CK-M subunits (i.e. 100% of CK-MM and 50% of CK-MB isozymes). The remaining activity, corresponding to CK-B fraction activity, is measured according to the IFCC reference method for measuring CK activity. CK-MB activity is then obtained by multiplying by 2 the remaining activity.

☛ **COMPOSITION**

Reagent 1 : R1

Imidazole buffer, pH 6.10	125 mmol/L
D-Glucose	25 mmol/L
N-Acetyl-L-Cysteine	25 mmol/L
Magnesium acetate	12.5 mmol/L
NADP	2.4 mmol/L
EDTA	2.0 mmol/L
Hexokinase	≥ 6 800 U/L
Sodium azide	< 0.1 %

The concentration of anti-CK-M antibody contained in reagent 1 is sufficient to inhibit 2000 U/L of CK-M at 37 °C.

Reagent 2 : R2

Imidazole buffer, pH 8.90	125 mmol/L
Creatine phosphate	250 mmol/L
ADP	15.2 mmol/L
AMP	23 mmol/L
Diadenosine pentaphosphate	103 µmol/L
G-6-PDH	≥ 8 800 U/L
Sodium azide	< 0.1 %

☛ **MATERIALS REQUIRED BUT NOT PROVIDED**

- CKMB-0900 CK-MB CONTROL
- General Laboratory equipment.
- Do not use materials that are not required as indicated above.

☛ **WARNINGS AND PRECAUTIONS**

- This *in vitro* diagnostic device is for professional use only.
- The reagents R1 and R2 are classified as hazardous (Imidazole):

DANGER : May damage the unborn child. Obtain special instructions before use. Wear protective gloves/ protective clothing/eye protection/ face protection. IF exposed or concerned: Get medical advice/attention.

- The reagents contain sodium azide which may react with lead or copper plumbing to form potentially explosive metal azides. When disposing of these reagents always flush with copious amounts of water to prevent azide buildup.
- For more information, refer to the Safety Data Sheet (SDS).
- Take normal precautions and adhere to good laboratory practice.
- Use clean or single use laboratory equipment only to avoid contamination.
- Do not interchange reagent vials from different kits.

☛ **STABILITIES**

Store at 2-8 °C and protect from light. Do not freeze.
Do not use after expiration dates indicated on the vial labels.

PREPARATION AND STABILITY OF WORKING REAGENT

- Preparation
 Mix in an empty vial (empty vial provided with the ref. CMSL-0230) 4 volumes of reagent R1 with 1 volume of reagent R2.

- Stability
 1 day at 20-25 °C
 2 weeks at 2-8 °C
On board stability :
 The on-board stability is specific for each analyzer. (Refer to § PERFORMANCE DATA).

➤PRODUCT DETERIORATION

- The product should be clear. Cloudiness would indicate deterioration.
 - Do not use the product if there is visible evidence of biological, chemical or physical deterioration.
 - Do not use the product if the damages of packaging might have an effect on the product performances (leakages, pierced vial).

SAMPLES (6)

Specimen
 - Serum free from hemolysis.
 - Do not use other specimens.
Warnings and precautions
 According to Good Laboratory Practice, sampling should be performed prior to the administration of drugs.
Storage and stability
 Samples must be analyzed immediately or stored protected from air and light 2 days at 2-8 °C or 1 month at -20 °C.

REFERENCE VALUES (3-5)

Serum (37 °C) : 0-25 U/L
 CK-MB activity must be compared to total CK activity (CK-MB/ total CK) x 100 < 6%

The following 3 factors are indicators of damage of cardiac muscle :
Total CK
 Men > 171 U/L
 Women > 145 U/L
CK-MB : > 25 U/L
Ratio : (CK-MB/ total CK) x 100: 6 - 25 %

Note : The quoted range should serve as a guide only. It is recommended that each laboratory verifies this range or establishes a reference interval for the intended population.

PROCEDURE

This reagent can be used on one-reagent procedure only
For ELITech Clinical Systems Selectra Analyzers, applications are available on request.
 Wavelength 340 nm
 Temperature: 37 °C
 Read against reagent blank.

Working Reagent	250 µL
Sample	10 µL

Mix and after a 130 second incubation, measure the change of absorbance per minute (ΔA/min.) during 159 seconds.

With Selectra TouchPro software, use the application included in the barcode available at the end of this insert.

CALCULATION

a) Total CK activity :
 Determination with CK NAC SL reagent.

b) CK-MB activity :
 At 340 nm, with a 1 cm light path cuvette :
 Activity (U/L) = ΔA/min. x 8 254

c) Percentage of CK-MB activity in sample :

$$\% \text{ CK-MB} = \frac{\text{CK-MB}}{\text{Total CK}} \times 100$$

Conversion factor : U/L x 0.0167 = µkat/L

QUALITY CONTROL

To check the accuracy of assays, control serum such as CK-MB CONTROL should be used. This control serum is prepared from unmodified human enzyme. The quality control must be performed and validated before the patient samples are assayed. The control frequency must be at least once a day, and should be adapted to Quality Control procedures of each laboratory and the regulatory requirements. Results should be within the defined ranges. If values fall outside of the defined ranges, each laboratory should take corrective measures. Quality control materials should be used in accordance with local guidelines.

WASTE MANAGEMENT

Disposal of all waste material should be in accordance with local, state and federal regulatory requirements.

PERFORMANCE DATA at 37 °C on ELITech Clinical Systems Selectra E Analyzers

- Measuring range
 The reagent is linear from 15 to 600 U/L.

- Limit of Detection (7)

Determined according to SFBC protocol, the detection limit is equal to 6 U/L.

- Precision

Within-run reproducibility	n	Mean	CV (%)
		U/L	
Normal level	20	11	16.3
Pathological level	20	152	1.0

Between-run reproducibility	n	Mean	CV (%)
		U/L	
Normal level	20	14	14.1
Pathological level	20	145	1.6

- Correlation

A comparative study was performed on this reagent between analyser ELITech Clinical Systems Selectra E and Cobas Mira on 60 human serum samples. The values cover the linearity range.
 The parameters of linear regression are as follows :
 Correlation coefficient: (r) = 0.9980
 Linear regression: y = 1.0007 x + 1.3 U/L

- Limitations/Interferences (2,7,8)

- Hemolysed sera should not be used since significant hemolysis may increase CK-MB concentration because of adenylate kinase release.

- The method will also measure any CK-BB isoenzyme present in serum alone or complexed with immunoglobulins (macro-CK). The activity of the isoenzyme is usually negligible, however, if a significant amount of CK-BB activity is present the CK-MB activity will be overestimated.

- Do not report results outside of the usable range.

- According to SFBC recommendations, some studies have been performed to determine the level of interference from different compounds :
Unconjugated bilirubin: Negative bias from 9 mg/dL (153.9 µmol/L) on human normal sera.
 Negative bias from 15 mg/dL (256.6 µmol/L) on human pathological sera.
Conjugated bilirubin: Negative bias from 2 mg/dL (34.2 µmol/L) on human normal sera.
 Negative bias from 6 mg/dL (102.6 µmol/L) on human pathological sera.
Turbidity: No significant interference up to 600 mg/dL (6.78 mmol/L) triglycerides equivalent.

- In very rare cases, monoclonal gammopathies (multiple myeloma), in particular IgM type (Waldenstrom's macroglobulinemia) can cause unreliable results.⁽⁹⁾
 - Many other substances and drugs may interfere. Some of them are listed in reviews published by Young.⁽¹⁰⁻¹¹⁾

- The results of this assay should only be interpreted in conjunction with other diagnostic test results, clinical findings and the patient's medical history.

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