CHOLESTEROL SL

► Références/References/ Composition du coffret/ Kit composition/
Referencias/ Referencias: Composición del kit/ Conteúdo da embalagem :

CHSL-0250	R	12 x	20 mL
CHSL-0455	R	6 x	45 mL
CHSL-0500	R	6 x	100 mL
CHSL-0700	R	4 x	250 mL
CHSL-0497	R	1 x	100 mL + Std 1 x 5 mL
CHSL-0507	R	6 x	100 mL + Std 1 x 5 mL
CHSL-0707	R	4 x	250 mL + Std 1 x 5 mL

►CONTROLE DE QUALIDADE

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BIBLIOGRAPHIE/BIBLIOGRAPHY

- Rifai, N., et al., *Lipids, Lipoproteins, and Apolipoproteins*. *Tietz Fundamentals of Clinical Chemistry*, 5th Ed., Burtis, C.A. & Ashwood, E.R. (W.B. Saunders eds. Philadelphia USA), (2001), 463.
- Naito, H.K., *Coronary Artery Disease and Disorders of Lipid Metabolism*. *Clinical Chemistry: Theory, Analysis, Correlation*, 4th Ed., Kaplan, L.A., Pesce, A.J., Kazmierczak, S.C. (Mosby, Inc. eds. St Louis USA), (2003), 603.
- Allain, C.C., et al., *Clin. Chem.*, (1974), 20, 470.
- Tietz, N.W. *Clinical guide to laboratory tests*, 3rd Ed., (W.B. Saunders eds. Philadelphia USA), (1995), 130.
- Expert Panel on Detection, *JAMA*, (2001), 285, 2486.
- Evaluation of the Linearity of the Measurement of Quantitative Procedures: a Statistical Approach; Approved Guideline. CLSI (NCCLS) document EP6-A (2003), 23 (16).
- Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline—Second Edition. CLSI (NCCLS) document EP5-A2 (2004), 24 (25).
- Method Comparison and Bias estimation Using Patient Samples; Approved Guideline—Second Edition. CLSI (NCCLS) document EP9-A2 (2002), 22 (19).
- Third Report of the NCEP Expert Panel on Detection, Evaluation, and Treatment of High Blood Cholesterol in Adults (Adult Treatment Panel III)-Final report. NIH Publication (2002), No. 02-5215.
- Interference Testing in Clinical Chemistry ; Approved Guideline—Second Edition. CLSI (NCCLS) document EP7-A2 (2005), 25 (27).
- Vassault, A., et al., *Ann. Biol. Clin.*, (1986), 44, 686.
- Berth, M. & Delanghe, J. *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature*. *Acta Clin Belg.*, (2004), 59, 263.
- Young, D.S., *Effects of preanalytical variables on clinical laboratory tests*, 2nd edition, AACC Press (1997).
- Young, D.S., *Effects of drugs on clinical laboratory tests*, 4th edition, AACC Press (1995).

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- Symbols used are defined on ISO-15223-1 standard, except those present below.

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CHOLESTEROL SL

WARNINGS AND PRECAUTIONS

- These *in vitro* diagnostic devices (reagent and standard) are for professional use only.
- The reagent R contains sodium azide which may react with lead or copper plumbing to form potentially explosive metal azides. When disposing of these reagents always flush with copious amounts of water to prevent azide buildup.
- Take normal precautions and adhere to good laboratory practice.
- Use clean or single use laboratory equipment only to avoid contamination.
- The standard should be immediately and tightly capped to prevent contamination and evaporation.
- For more information, Safety Data Sheet (SDS) is available on request for professional user.

STABILITIES

Store at 2-8 °C and protect from light. Do not freeze. Do not use after expiration dates indicated on the vial labels.

On board stability :

The on-board stability is specific for each analyzer. (Refer to § PERFORMANCE DATA).

PREPARATION

The reagent and standard are ready to use.

PRODUCT DETERIORATION

- The product should be clear. Cloudiness would indicate deterioration.
- Do not use the product if there is visible evidence of biological, chemical or physical deterioration.
- Do not use the product if the damages of packaging might have an effect on the product performances (leakages, pierced vial)

SAMPLES (1,4)

Specimen

- Serum or lithium heparinized plasma from fasting patients (≥ 12 hours).

- Do not use other specimens.

Warnings and precautions

- According to Good Laboratory Practice, sampling should be performed prior to the administration of drugs. Sampling could lead to false results if performed during or immediately after the administration of some drugs.
- Separate from cells within 2 hours.
- To reduce biological variability, collection of samples should follow standardized conditions as recommended by NCEP.

Storage and stability

- Samples are stable 5 to 7 days if stored at 2-8 °C, 3 months at -20 °C and several years at -70 °C. Avoid repeated freezing and thawing.

REFERENCE VALUES (5)

The NCEP (American National Cholesterol Education Program) has established the following classification for total cholesterol levels according to the risk of developing coronary heart disease:

Risk Classification

	Level (mg/dL)	Level (mmol/L)
Normal	< 200	< 5.17
Borderline high	200 - 239	5.17 - 6.18
High	≥ 240	≥ 6.21

Note : The quoted range should serve as a guide only. It is recommended that each laboratory verifies this range or establishes a reference interval for the intended population.

PROCEDURE

Manual Procedure

Wavelength : 505 nm
Optical path : 1 cm
Sample/Reagent ratio : 1:100
Temperature : 37 °C
Read against reagent blank.

	CALIBRATION	TEST
Reagent R	1 000 µL	1 000 µL
Standard/Calibrator	10 µL	-
Sample	-	10 µL

Mix and read the absorbances (A) after an incubation of 5 minutes.

Automatic Procedure

These reagents may be used in several automatic analyzers. For ELITech Selectra Analyzers, validated applications are available on request. For Selectra TouchPro software, use the application included in the barcode available at the end of this insert.

CALCULATION

A Sample $\times n$ n = Calibrator/ standard concentration

Conversion factor : mg/dL x 0.0259 = mmol/L mg/dL x 0.01 = g/L

CALIBRATION

For the references CHSL-0497/0507/0707: For calibration, use either multiparametric calibrator ELICAL 2 or Cholesterol Standard 200 mg/dL.
For the references CHSL-0250/0455/0500/0700: For calibration, use multiparametric calibrator ELICAL 2.

Concentration values of Cholesterol Standard 200 mg/dL and multiparametric calibrator ELICAL 2 are traceable to ID-GC-MS (Isotope Dilution - Gas Chromatography - Mass Spectrometry) reference method.

Calibration frequency : The calibration is specific for each analyzer. (Refer to § PERFORMANCE DATA).

QUALITY CONTROL

To check the accuracy of assays, control sera such as ELITROL I and ELITROL II should be used. These controls must be performed and validated before the patient samples are assayed. The control frequency must be at least once a day, after each calibration and should be adapted to Quality Control procedures of each laboratory and the regulatory requirements. Results should be within the defined ranges. If values fall outside of the defined ranges, each laboratory should take corrective measures. Quality control materials should be used in accordance with local guidelines.

WASTE MANAGEMENT

Disposal of all waste material should be in accordance with local, state and federal regulatory requirements.

PERFORMANCE DATA at 37 °C on ELITech Clinical Systems Selectra ProM Analyzers

- Measuring range

Determined according to CLSI⁽⁶⁾ EP6-A protocol, the measuring range is from 20 to 600 mg/dL (0.52 to 15.52 mmol/L). Samples having greater concentrations should be diluted 1:5 with NaCl 9 g/L solution and re-assayed. This procedure extends the measuring range up to 3000 mg/dL (77.59 mmol/L).

For users with Selectra TouchPro software, the «dilute» function performs the sample dilution automatically. Results take the dilution into account.

Precision

Determined according to CLSI EP5-A2 protocol⁽⁷⁾.

	Mean	Within-run	Total
n	mg/dL	mmol/L	CV (%)
Low level	80	115	2.97
Medium level	80	184	4.76
High level	80	292	7.55

Correlation

A comparative study has been performed between an ELITech Clinical Systems Selectra ProM Analyzer and another FDA-approved system equipment (enzymatic and colorimetric method) on 100 human serum samples according to CLSI EP9-A2 protocol⁽⁸⁾.

The values were between 20 and 575 mg/dL (0.52 and 14.87 mmol/L). The parameters of the linear regressions are as follows :

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Linear regression: $y = 1.016 x + 0 \text{ mg/dL}$

Limitations/Interferences

- At least two measurements of cholesterol should be made on separate occasions before any medical decision.⁽⁹⁾

- Do not report results outside of the usable range.

- Studies have been performed to determine the level of interference from different compounds according CLSI EP7-A2 protocol⁽¹⁰⁾ and SFBC recommendations⁽¹¹⁾. Recovery is within $\pm 10\%$ of initial value of cholesterol concentration of 116 mg/dL and 309 mg/dL.

COMPOSIÇÃO

Reactivos : R

Tampón de Good, pH 6,7

Fenol 24 mmol/L

Amino-4-antipirina 0,5 mmol/L

Colesterol esterase ≥ 180 U/L

Colesterol oxidase ≥ 200 U/L

Peroxidasa $\geq 1\ 000$ U/L

Azida sódica < 0,1 %

Media 1 mg/dL

Intra-série 1,1

Total 2,1

CV (%) 1,1

Unconjugated Bilirubin: No significant interference up to 6.0 mg/dL (103 µmol/L).

Conjugated Bilirubin: No significant interference up to 5.9 mg/dL (101 µmol/L).

Hemoglobin : No significant interference up to 300 mg/dL.

Turbo: No significant interference up to 614 mg/dL (6.94 mmol/L) triglyceride equivalent.

Ascorbic acid : No significant interference up to 4.0 mg/dL. Concentraciones above the therapeutic levels will interfere and cause erroneous results.

Methyl-dopa: No significant interference up to 1.6 mg/dL. Concentraciones above the therapeutic levels will interfere and cause erroneous results.

Uric acid: No significant interference up to 23.7 mg/dL (1410 µmol/L).

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