## AMOSTRAS <sup>(1,3,4)</sup>

### Amostras

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2. Daly J.A., Ertingshausen G., *Direct method for determining inorganic phosphorus in serum with the Centrifchem. Clin. Chem.*, (1972), **18**, 263.

3. Itani, O., Tsang, R.C., *Bone disease. Clinical Chemistry: Theory Analysis. Correlation*, 5<sup>th</sup> Ed., Kaplan, L.A. Pesce, (Mosby Inc. eds St Louis USA), (2010), 614 and appendix.

4. Wu, A. H. B., *Clinical guide to laboratory tests*, 4<sup>th</sup> Ed., (W.B. Saunders eds. Philadelphia USA), (2006), 852.

5. *Evaluation of the Linearity of the Measurement of Quantitative Procedures: a Statistical Approach; Approved Guideline.* CLSI (NCCLS) document EP6-A (2003), **23** (16).

6. *Protocols for Determination of Limits of Detection and Limits of Quantification; Approved Guideline.* CLSI (NCCLS) document EP17-A (2004), **24** (34).

7. *Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline - Second Edition.* CLSI (NCCLS) document EP5-A2 (2004), **24** (25).

8. *Method Comparison and Bias estimation Using Patient Samples; Approved Guideline - Second Edition.* CLSI (NCCLS) document EP9-A2 (2002), **22** (19).

9. *Interference Testing in Clinical Chemistry ; Approved Guideline - Second Edition.* CLSI (NCCLS) document EP7-A2 (2005), **25** (27).

10. Vassault A., et al., *Ann. Biol. Clin.*, (1986), **44**, 686.

11. Berth, M. & Delanghe, J., *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature.* *Acta Clin. Belg.*, (2004), **59**, 263.

12. Young, D. S., *Effects of preanalytical variables on clinical laboratory tests*, 2<sup>nd</sup> Ed., AACCC Press, (1997).

13. Young, D. S., *Effects of drugs on clinical laboratory tests*, 4<sup>th</sup> Ed., AACCC Press, (1995).

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Symbols used are defined on ISO-15223-1 standard, except those presented below.

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**PRINCIPLE** <sup>(2)</sup>

Determination of inorganic phosphorus is made according to the following reaction:

**COMPOSITION**

<b>Reagent: R</b>	
Sulfuric acid	210 mmol/L
Ammonium molybdate	650 µmol/L
<b>Standard: Std</b>	
Phosphorus	5 mg/dL 1.61 mmol/L

**MATERIALS REQUIRED BUT NOT PROVIDED**

- CONT-0060 ELITROL I
- CONT-0160 ELITROL II
- Normal saline solution (NaCl 9 g/L).
- General Laboratory equipment.
- Do not use materials that are not required as indicated above.

**☛ WARNINGS AND PRECAUTIONS**

- These *in vitro* diagnostic devices (reagent and standard) are for professional use only.
- The reagent R is classified as hazardous :

**WARNING :** Causes skin irritation. Causes serious eye irritation. Wear protective gloves / eye protection / face protection. IF IN EYES : Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists: Get medical advice/attention.

- For more information, refer to the Safety Data Sheet (SDS).

- The standard should be immediately and tightly capped to prevent contamination and evaporation.
- Take normal precautions and adhere to good laboratory practice.

- Use clean or single use laboratory equipment only to avoid contamination.

**STABILITIES**

**Store at 2-25 °C and protect from light. Do not freeze.** Do not use after expiration dates indicated on the vial labels.

**On board stability :**

The on-board stability is specific for each analyzer. (Refer to § PERFORMANCE DATA).

**PREPARATION**

The reagent and standard are ready to use.

**PRODUCT DETERIORATION**

- The product should be clear. Cloudiness would indicate deterioration.
- Do not use the product if there is visible evidence of biological, chemical or physical deterioration.
- Do not use the product if the damages of packaging might have an effect on the product performances (leakages, pierced vial).

**SAMPLES** <sup>(1,3,4)</sup>

**Specimen**

- Serum free of hemolysis from fasting patient.
  - Lithium heparinized plasma free of hemolysis from fasting patient.
  - Urine to acidify (pH<3) after collection and to dilute at 1/10 with saline solution NaCl 9 g/L before analysis (when there is no predilution by the analyser). For users with Selectra Touch Pro software, Pre-dilution is performed automatically. Results take the dilution into account.
  - Do not use other specimens.
- Warnings and precautions**
- According to Good Laboratory Practice, sampling should be performed prior to the administration of drugs.
  - Separate from erythrocytes promptly after collection.
  - Do not use visibly turbid or icteric samples.

**Storage**

- Plasma and serum are stable 4 days at room temperature, 7 days at 2-8 °C, and several months at -20 °C.
- Acidified urines are stable for 6 months.

**REFERENCE VALUES** <sup>(1,4)</sup>

Serum, Plasma :	2.5 - 4.5 mg/dL 0.81 - 1.45 mmol/L
Urine:	0.4 - 1.3 g/24 h 12.9 - 42.0 mmol/24 h 26.7 - 86.7 mg/dL* 8.6 - 28.0 mmol/L*

\* for an urinary volume of 1.5 L per 24 hours

**Note :** The quoted range should serve as a guide only. It is recommended that each laboratory verifies this range or establishes a reference interval for the intended population.

**☛PROCEDURE**

**Manual Procedure**

**Wavelength :** 340 nm  
**Optical path :** 1 cm  
**Sample/reagent ratio :** 1:100  
**Temperature :** 37 °C  
*Read against reagent blank.*

	CALIBRATION	TEST
<b>Reagent R</b>	1 000 µL	1 000 µL
<b>Standard</b>	10 µL	-
<b>Sample</b>	-	10 µL

Mix and read the absorbances (A) after an incubation of 5 minutes.

**Automatic Procedure**

These reagents may be used in several automatic analyzers. For ELITech Selectra Analyzers, validated applications are available on request. For Selectra TouchPro software, use the application included in the barcode available at the end of this insert.

**CALCULATION**

$$\frac{(A) \text{ Sample}}{(A) \text{ Standard}} \times n \quad n = \text{standard concentration}$$

Take dilution factor into account for the calculation of phosphorus concentration in urine.

Conversion factor : mg/dL x 0.3229 = mmol/L

**CALIBRATION**

Phosphorus standard is traceable to the Standard Reference Material NIST SRM 3139a.

**Calibration frequency :** The calibration is specific for each analyzer. (Refer to § PERFORMANCE DATA).

**QUALITY CONTROL**

To check the accuracy of assays, control sera such as ELITROL I and ELITROL II should be used. These controls must be performed and validated before the patient samples are assayed. The control frequency must be at least once a day, after each calibration and should be adapted to Quality Control procedures of each laboratory and the regulatory requirements. Results should be within the defined ranges. If values fall outside of the defined ranges, each laboratory should take corrective measures. Quality control materials should be used in accordance with local guidelines.

**WASTE MANAGEMENT**

Disposal of all waste material should be in accordance with local, state and federal regulatory requirements.

**PERFORMANCE DATA at 37 °C on ELITech Clinical Systems Selectra ProM Analyzers**

**- Measuring range**  
Determined according to CLSI EP6-A protocol<sup>(5)</sup>.

**a) Serum/Plasma**  
The measuring range is from 2.0 to 20.0 mg/dL (0.65 to 6.46 mmol/L).

**☛ b) Urine**  
The measuring range is from 10.0 to 200.0 mg/dL (3.2 to 64.6 mmol/L). Samples having greater concentrations should be diluted 1:5 with NaCl 9 g/L solution and re-assayed. This procedure extends the measuring range up to 1000.0 mg/dL (322.9 mmol/L).

For Selectra TouchPro users, the «dilute» function performs the sample dilution automatically. Results take the dilution into account.

**- Limit of Detection (LoD) and Limit of Quantification (LoQ)**  
Determined according to CLSI EP17-A protocol<sup>(6)</sup>.

**a) Serum/Plasma**  
LoD = 0.12 mg/dL (0.04 mmol/L)  
LoQ = 1.00 mg/dL (0.32 mmol/L)

**b) Urine**  
LoD = 0.2 mg/dL (0.1 mmol/L)  
LoQ = 10.0 mg/dL (3.2 mmol/L)

**- Precision**  
Determined according to CLSI EP5-A2 protocol<sup>(7)</sup>.

**a) Serum/Plasma**

		Mean	Within-run	Total	
	n	mg/dL	mmol/L	CV (%)	
<b>Low level</b>	80	2.39	0.77	0.6	2.1
<b>Medium level</b>	80	4.77	1.54	0.4	2.3
<b>High level</b>	80	9.02	2.91	0.7	2.3

**b) Urine**

		Mean	Within-run	Total	
	n	mg/dL	mmol/L	CV (%)	
<b>Low level</b>	80	25.2	8.1	1.0	2.0
<b>Medium level</b>	80	72.5	23.4	0.6	1.0
<b>High level</b>	80	144.2	46.6	0.5	1.2

**- Correlation**

**a) Serum/Plasma**

A comparative study has been performed between an ELITech Clinical Systems Selectra ProM Analyzer and another FDA-approved system equipment (molybdate-UV method) on 100 human sera samples according to CLSI EP9-A2 protocol<sup>(8)</sup>. The values were between 2.01 and 18.82 mg/dL (0.65 and 6.08 mmol/L). The parameters of the linear regressions are as follows:  
Correlation coefficient : (r) = 0.999  
Linear regression: y = 1.004 x + 0.05 mg/dL (0.02 mmol/L)

**b) Urine**

A comparative study has been performed between an ELITech Clinical Systems Selectra ProM Analyzer and another FDA-approved system equipment (molybdate-UV method) on 52 urine samples according to CLSI EP9-A2 protocol<sup>(8)</sup>. The values were between 10.2 and 182.3 mg/dL (3.3 and 58.9 mmol/L). The parameters of the linear regressions are as follows:  
Correlation coefficient : (r) = 0.999  
Linear regression: y = 0.883 x - 0.3 mg/dL (0.1 mmol/L)

**- Limitations/interferences**

- Do not report results outside of the usable range.

**a) Serum plasma**

Studies have been performed to determine the level of interference from different compounds according to CLSI EP7-A2 protocol<sup>(9)</sup> and SFBC recommendations<sup>(10)</sup>. Recovery is within ±10% of initial value at phosphorus concentration of 2.17 mg/dL and 4.65 mg/dL.

**Unconjugated bilirubin:** No significant interference up to 15 mg/dL (257 µmol/L).  
**Conjugated bilirubin:** No significant interference up to 3.0 mg/dL (50 µmol/L).  
**Glucose :** No significant interference up to 528 mg/dL (29.31 mmol/L).  
**Triglycerides:** No significant interference up to 729 mg/dL (8.24 mmol/L).  
**Turbidity:** Interference occurs at all levels of Intralipid®.  
**Hemoglobin:** No significant interference up to 50 mg/dL.  
**Ascorbic acid:** No significant interference up to 20.0 mg/dL.  
**Methyl-Dopa:** No significant interference up to 1 mg/dL.

**- In very rare cases, monoclonal gammopathies (multiple myeloma), in particular IgM type (Waldenström's macroglobulinemia) can cause unreliable results.**<sup>(11)</sup>

- Many other substances and drugs may interfere. Some of them are listed in reviews published by Young.<sup>(12-13)</sup>

- The results of this assay should only be interpreted in conjunction with other diagnostic test results, clinical findings and the patient's medical history.

**b) Urine**  
Studies have been performed to determine the level of interference from different compounds according to CLSI EP7-A2 protocol<sup>(9)</sup>. Recovery is within ±10% of initial value at phosphorus concentration of 25 mg/dL and 150 mg/dL.

**Conjugated bilirubin:** No significant interference up to 22.1 mg/dL (378 µmol/L).  
**Hemoglobin:** No significant interference up to 500 mg/dL.  
**Uric acid :** No significant interference up to 120 mg/dL (38.7 mmol/L).

**Urea :** No significant interference up to 6000 mg/dL (999 mmol/L).  
**Ascorbic acid:** No significant interference up to 20 mg/dL.

- Many other substances and drugs may interfere. Some of them are listed in reviews published by Young.<sup>(12-13)</sup>

- The results of this assay should only be interpreted in conjunction with other diagnostic test results, clinical findings and the patient's medical history.

**- On board stability/Calibration frequency**  
**On Board Stability:** 28 days  
**Calibration frequency:** 28 days  
Recalibrate when reagent lots change, when quality control results fall outside the established range, and after a maintenance operation.

**☛ These performances have been obtained using ELITech Selectra ProM analyzer. Results may vary if a different instrument or a manual procedure is used. The performances of applications not validated by ELITech are not warranted and must be defined by the user.**

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