➤ **Références/Referenes/ Referencias/ Referências:**

GPSL-0497
GPSL-0507
GPSL-0707
GPSL-0250
GPSL-0455
GPSL-0500
GPSL-0700

Composition du coffret/ Kit composition/ Composición del kit/ Conteúdo da embalagem :

R 1 x 100 mL + Std 1 x 5 mL
R 6 x 100 mL + Std 1 x 5 mL
R 4 x 250 mL + Std 1 x 5 mL
R 12 x 20 mL
R 6 x 45 mL
R 6 x 100 mL
R 4 x 250 mL

➤ **SYMBOLS/SYMBOLS/ SÍMBOLOS/SÍMBOLOS**

- Les symboles utilisés sont décrits dans la norme ISO-15223-1 hormis ceux présentés ci-dessous.
- Symbols used are defined on ISO-15223-1 standard, except those presented below.
- Los símbolos utilizados son descritos en la norma ISO-15223-1 a la excepción de los presentados a continuación.
- Os símbolos utilizados são definidos na norma ISO-15223-1, exceto os apresentados abaixo.

CONT	Contient Content Contiene Conteúdo
R	Réactif Reagent Reactivo Reagente
Std	Standard Standard Estándar Padrão
CE	Conformité Européenne European Conformity Conformidad Europea Conformidade Europeia

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BIBLIOGRAPHIE/BIBLIOGRAPHY BIBLIOGRAFIA/BIBLIOGRAFIA

- Sacks, D.B., *Carbohydrates. Tietz Fundamentals of Clinical Chemistry*, 6th Ed., Burtis, C.A. & Ashwood, E.R. (W.B. Saunders eds. Philadelphia USA), (2008), 373.
- Dods, R.F., *Diabetes Mellitus. Clinical Chemistry : Theory, Analysis, Correlation*, 5th Ed., Kaplan, L.A, Pesce, A.J., S.C., (Mosby Inc. eds St Louis USA), (2010), 729 and appendix
- Wu, H.B., *General Clinical Tests. Tietz Clinical guide to laboratory tests*, 4th Ed., (W.B. Saunders eds. Philadelphia USA), (2006), 444.
- Trinder, P., *Determination of glucose in blood using glucose oxidase with an alternative oxygen acceptor. Ann. Clin. Biochem.*, (1969), **6**, 24.
- Guder W.G, *Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma, and serum samples, World Health Organization , WHO/ DIL/LAB/99.1 Rev.2, (2002).*
- Evaluation of the Linearity of the Measurement of Quantitative Procedures: a Statistical Approach; Approved Guideline.* CLSI(NCCLS) document EP6-A (2003), **23** (16).
- Protocols for Determination of Limits of Detection and Limits of Quantification; Approved Guideline.* CLSI (NCCLS) document EP17-A (2004), **24** (34).
- Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline - Second Edition.* CLSI (NCCLS) document EP5-A2 (2004), **24** (25).
- Method Comparison and Bias estimation Using Patient Samples; Approved Guideline - Second Edition.* CLSI (NCCLS) document EP9-A2 (2002), **22** (19).
- Interference Testing in Clinical Chemistry ; Approved Guideline - Second Edition.* CLSI (NCCLS) document EP7-A2 (2005), **25** (27).
- Berth, M. & Delanghe, J. *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature, Acta Clin Belg.*, (2004), **59**, 263.
- Young, D. S., *Effects of preanalytical variables on clinical laboratory tests*, 2nd Ed., AACCC Press, (1997).
- Young, D. S., *Effects of drugs on clinical laboratory tests*, 4th Ed., AACCC Press, (1995).

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- Only for ref. **GPSL-0250/0455**, used with Selectra TouchPro software.

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English - EN

➤ **INTENDED USE**

ELITech Clinical Systems GLUCOSE PAP SL is an *in vitro* diagnostic reagent intended for the quantitative determination of glucose in human serum and plasma samples.

CLINICAL SIGNIFICANCE (1-3)


Glucose is the main source of energy for the human body. Glucose is converted either into glycogen to be stocked in the liver or into triglycerides to be stocked in fatty tissues. Glucose concentration in blood is regulated by several hormones, including two antagonists : insulin and glucagon.
Quantification of glucose in blood is used to diagnose metabolic carbohydrates disorders such as diabetes, idiopathic hypoglycemia and pancreatic disease. The main physiological troubles are linked to hyperglycaemia (type I Diabetes mellitus and type II Diabetes mellitus). Type I diabetes mellitus is insulin-dependent, and appears mainly before 30 years old. Type II diabetes mellitus is non-insulin-dependent, and usually appears after 40 years old, but can occur earlier for obese people. Other diabetes have secondary origin, and appear after endocrinal or hepatic diseases.

METHOD (4)

Enzymatic colorimetric.
Trinder - End Point.

➤ **PRINCIPE (4)**

Enzymatic determination of glucose according to the following reactions:

❖ Références/Referenes/ Referencias/ Referências:

GPSL-0497
GPSL-0507
GPSL-0707
GPSL-0250
GPSL-0455
GPSL-0500
GPSL-0700

Composition du coffret/ Kit composition/ Composición del kit/ Conteúdo da embalagem :

R 1 x 100 mL + **Std** 1 x 5 mL
R 6 x 100 mL + **Std** 1 x 5 mL
R 4 x 250 mL + **Std** 1 x 5 mL
R 12 x 20 mL
R 6 x 45 mL
R 6 x 100 mL
R 4 x 250 mL

- Use clean or single use laboratory equipment only to avoid contamination.

- The standard should be immediately and tightly capped to prevent contamination and evaporation.

- For more information, Safety Data Sheet (SDS) is available on request for professional user.

❖ STABILITIES

Store at 2-8 °C and protect from light. Do not freeze. Do not use after expiration dates indicated on the vial labels.

On board stability :

The on-board stability is specific for each analyzer. (Refer to § PERFORMANCE DATA).

❖ PREPARATION

The reagent and standard are ready to use.

❖ PRODUCT DETERIORATION

- The product should be clear. Cloudiness would indicate deterioration.

- Do not use the product if there is visible evidence of biological, chemical or physical deterioration.

- Do not use the product if the damages of packaging might have an effect on the product performances (leakages, pierced vial).

SAMPLES ^(1,5)

Specimen

- Serum.

- Plasma in sodium fluoride/potassium oxalate (glycolysis inhibitor) or lithium heparin. Rather, it is recommended to use the plasma collected on sodium fluoride/potassium oxalate.

- Samples must be free from haemolysis.

- Do not use other specimens.

❖ Warnings and precautions

- According to Good Laboratory Practice, sampling should be performed prior to the administration of drugs. Sampling could lead to false results if performed during or immediately after the administration of some drugs.

- Samples collected without sodium fluoride must be separated from cells and clot promptly after collection to minimize loss of glucose through glycolysis (decrease of 5-7% in one hour in whole blood at room temperature).

Storage and stability

- Serum and plasma collected without sodium fluoride are stable for 8 hours at room temperature and up to 3 days at 2-8 °C.

- Plasmas collected with sodium fluoride are stable for 2 days at room temperature and up to 7 days at 2-8 °C.

REFERENCE VALUES ⁽⁸⁾

Serum, plasma : 74 - 106 mg/dL
4.1 - 5.9 mmol/L

Note : The quoted range should serve as a guide only. It is recommended that each laboratory verifies this range or establishes a reference interval for the intended population.

❖ PROCEDURE

Manual Procedure

Wavelength : 505 nm
Optical path : 1 cm
Sample/reagent ratio : 1:100
Temperature: 37 °C
Read against reagent blank.

	CALIBRATION	TEST
Reagent R	1 000 µL	1 000 µL
Standard/ Calibrator	10 µL	-
Sample	-	10 µL

Mix and read the absorbances (A) after an incubation of 10 minutes.

Automatic Procedure

These reagents may be used in several automatic analyzers. For ELITech Selectra Analyzers, validated applications are available on request. For Selectra TouchPro software, use the application included in the barcode available at the end of this insert.

❖ CALCULATION

$A_{\text{Sample}} \times n \quad n = \text{Calibrator/ standard concentration}$
 $A_{\text{Calibrator/ Standard}}$

Conversion factor : mg/dL x 0.0555 = mmol/L

❖ CALIBRATION

For the references GPSL-0497/0507/0707: For calibration, use either multiparametric calibrator ELICAL 2 or Glucose Standard 100 mg/dL.

For the references GPSL-0250/0455/0500/0700: For calibration, use multiparametric calibrator ELICAL 2.

Concentration values of Glucose Standard 100 mg/dL and multiparametric calibrator ELICAL 2 are traceable to the reference method ID-MS (isotope Dilution - Mass Spectrometry)

Calibration frequency : The calibration is specific for each analyzer. (Refer to § PERFORMANCE DATA).

❖ QUALITY CONTROL

To check the accuracy of assays, control sera such as ELITROL I and ELITROL II should be used. These controls must be performed and validated before the patient samples are assayed. The control frequency must be at least once a day, after each calibration and should be adapted to Quality Control procedures of each laboratory and the regulatory requirements. Results should be within the defined ranges. If values fall outside of the defined ranges, each laboratory should take corrective measures. Quality control materials should be used in accordance with local guidelines.

WASTE MANAGEMENT

Disposal of all waste material should be in accordance with local, state and federal regulatory requirements.

PERFORMANCE DATA at 37 °C on ELITech Clinical Systems Selectra ProM Analyzers

❖ Measuring range

Determined according to CLSI⁽⁹⁾ EP6-A protocol, the measuring range is from 20.0 to 400.0 mg/dL (1.11 to 22.20 mmol/L).

Samples having greater concentrations should be diluted 1:5 with NaCl 9 g/L solution and re-assayed. This procedure extends the measuring range up to 2000.0 mg/dL (111.01 mmol/L).

For users with Selectra TouchPro software, the «dilute» function performs the sample dilution automatically. Results take the dilution into account.

Limit of Detection (LoD) and Limit of Quantification (LoQ)

Determined according to CLSI EP17-A protocol⁽⁷⁾.
LoD = 0.2 mg/dL (0.01 mmol/L)
LoQ = 10.0 mg/dL (0.56 mmol/L)

- Precision

Determined according to CLSI EP5-A2 protocol⁽⁸⁾.

		Mean		Within-run	Total
	n	mg/dL	mmol/L	CV (%)	
Low level	80	37.4	2.08	0.7	1.6
Medium level	80	113.1	6.28	0.5	0.9
High level	80	284.0	15.76	0.7	1.3

- Correlation

A comparative study has been performed between an ELITech Clinical Systems Selectra ProM Analyzer and another FDA-approved system equipment (Glucose oxidase method) on 100 human serum samples according to CLSI EP9-A2 protocol⁽⁹⁾. The values were between 22.2 and 384.9 mg/dL (1.23 and 21.36 mmol/L). The parameters of the linear regressions are as follows :

Correlation coefficient: (r) = 1.000
Linear regression: y = 0.989 x + 1.1 mg/dL (0.06 mmol/L)

- Limitations/Interferences

- Do not report results outside of the usable range.

- Studies have been performed to determine the level of interference from different compounds according CLSI EP7-A2 protocol⁽¹⁰⁾. Recovery is within ±10% of initial value of glucose concentration of 36.0 mg/dL, 108.1 mg/dL and 400.0 mg/dL.
Unconjugated bilirubin: No significant interference up to 6.0 mg/dL (103 µmol/L).

Conjugated bilirubin: No significant interference up to 5.9 mg/dL (101 µmol/L).

Hemoglobin: No significant interference up to 300 mg/dL.

Triglycerides: No significant interference up to 920 mg/dL (10.40 mmol/L).

Ascorbic acid: No significant interference up to 2.0 mg/dL.

Uric acid: No significant interference up to 23.0 mg/dL (1368 µmol/L).

Methyl dopa: No significant interference up to 0.8 mg/dL.

L-Dopa: Induces falsely low results at therapeutic concentrations.

Tolazamide: No significant interference up to 40.0 mg/dL.

Acetaminophen: No significant interference up to 30.0 mg/dL.

- In very rare cases, monoclonal gammopathies (multiple myeloma), in particular IgM type (Waldenström's macroglobulinemia) can cause unreliable results.⁽¹¹⁾

- Results can be falsely lowered by significant levels in sample of NAC (N-Acetyl-Cysteine), NAPQI (metabolite of acetaminophene (paracetamol)) or metemazole.

- Many other substances and drugs may interfere. Some of them are listed in reviews published by Young.^(12,13)

- The results of this assay should only be interpreted in conjunction with other diagnostic test results, clinical findings and the patient's medical history.

- On board stability/Calibration frequency

On Board Stability: 28 days
Calibration frequency: 28 days
Recalibrate when reagent lots change, when quality control results fall outside the established range, and after a maintenance operation.

❖ These performances have been obtained using ELITech Selectra ProM analyzer. Results may vary if a different instrument or a manual procedure is used. The performances of applications not validated by ELITech are not warranted and must be defined by the user.

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