HbA1c

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BIBLIOGRAPHIE/BIBLIOGRAPHY

BIBLIOGRAFIA/BIBLIOGRAFIA

1. Sacks, D.B., Carbohydrates, Tietz Fundamentals of Clinical Chemistry, 5th Ed., Burtis, C.A. & Ashwood, E.R. (W.B. Saunders eds. Philadelphia USA), (2001), 427.

2. Sacks, D.B., et al., Clinical Chemistry, (2002), 48, 436.

3. Guder, W.G., et al., Use of Anticoagulants in Diagnostic Laboratory Investigations and stability of blood, plasma and serum samples, WHO/DIL/LAB/99.1 Rev.2. 2002.

4. American Diabetes Association, Standards of Medical Diabetes Care in Diabetes-2008, Diabetes Care, (2008), 31 Suppl 1, S12.

5. Weykamp, C., et al., A review of the Challenge in Measuring Hemoglobin A1c, J Diabetes Sci Technol, (2009), 3, 439.

6. Evaluation of the Linearity of the Measurement of Quantitative Procedures: a Statistical Approach; Approved Guideline. CLSI (NCCLS) document EP6-A (2003), 23 (16).

7. Protocols for Determination of Limits of Detection and Limits of Quantification; Approved Guideline. CLSI (NCCLS) document EP17-A (2004), 24 (34).

8. Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline—Second Edition. CLSI (NCCLS) document EP5-A2 (2004), 24 (25).

9. Method Comparison and Bias estimation Using Patient Samples; Approved Guideline—Second Edition. CLSI (NCCLS) document EP9-A2 (2002), 212 (19).

10. Goldstein DE, Little RR, Wiedmeyer HM, England JD, and McKenzie EM. Glycated hemoglobin: methodologies and clinical applications. Clin. Chem. (1986) 32, B64-B70.

11. Interference Testing in Clinical Chemistry ; Approved Guideline - Second Edition. CLSI (NCCLS) document EP7-A2 (2005), 25 (27).

12. Vassault, A., et al., Ann. Biol. Clin., (1986), 44, 686.

13. Young, D.S., Effects of preanalytical variables on clinical laboratory tests, 2nd Ed., AACC Press, (1997).

14. Young, D.S., Effects of drugs on clinical laboratory tests, 4th Ed., AACC Press, (1995).

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PRINCIPLE

1st reaction:

The sample is mixed with R1 which contains uncoated latex particles. As total haemoglobin and HbA1c have the same absorption affinity for these particles, the %HbA1c present in the sample is proportional to latex-bound HbA1c.

2nd reaction:

Reagent R2 contains a mouse anti-human HbA1c monoclonal antibody and a goat anti-mouse polyclonal antibody. Agglutination complexes are formed from the interaction between latex-bound HbA1c and the corresponding antibodies. Turbidity created by these aggregates is proportional to the amount of latex-bound HbA1c and therefore is proportional to the %HbA1c in the sample. A non-linear calibration curve is used to obtain the %HbA1c.

COMPOSITION

Reagent 1: R1

Suspended latex particles 0.13 %
Sodium azide < 0.1 %

Buffer, stabilizers

Working reagent R2: R2a+R2b
Mouse anti-human HbA1c monoclonal antibody 0.05 mg/mL

Goat anti-mouse IgG polyclonal antibody 0.08 mg/dL
Sodium azide < 0.1 %

Buffer, stabilizers

Reagent 3: R3 (hemolysis reagent)
Sodium azide < 0.1 %

Water, stabilizer

MATERIALS REQUIRED BUT NOT PROVIDED

- HBAC-0043 HbA1c CALIBRATOR SET
- HBAC-0049 HbA1c CONTROL L & H
- Normal saline solution (NaCl 9 g/L).
- General Laboratory equipment.

- Do not use materials that are not required as indicated above.

WARNINGS AND PRECAUTIONS

- This *in vitro* diagnostic device is for professional use only.
- The reagents R1, R2a, R2b et R3 contain sodium azide which may react with lead or copper plumbing to form potentially explosive metal azides. When disposing of these reagents always flush with copious amounts of water to prevent azide buildup.
- Take normal precautions and adhere to good laboratory practice.
- Use clean or single use laboratory equipment only to avoid contamination.
- Do not interchange reagent vials from different kits.
- For more information, Safety Data Sheet (SDS) is available on request for professional user.

STABILITIES

Store at 2-8 °C and protect from light. Do not freeze.
Do not use after expiration dates indicated on the vial labels.

On board stability :
The on-board stability is specific for each analyzer.
(Refer to § PERFORMANCE DATA).

PREPARATION

Reagent R1: Ready for use
Working reagent R2: Pour entire content of R2b in R2a. Mix gently.

Stability: 4 weeks at 2-8 °C.
Reagent R3: Ready for use.

PRODUCT DETERIORATION

- The solutions of reagents R2a, R2b and R3 should be clear. Cloudiness would indicate deterioration. The reagent R1 solution may present slightly hazy appearance without consequence on the performances of the product.
- Do not use the product if there is visible evidence of biological, chemical or physical deterioration.
- Do not use the product if the damages of packaging might have an effect on the product performance (leakages, pierced vial).

SAMPLES (2,3)

Specimen
- Whole blood collected on EDTA.
- Do not use other specimens.

Warnings and precautions
According to Good Laboratory Practice, sampling should be performed prior to the administration of drugs.

Storage and stability

Whole blood samples are stable 3 days at 20-25 °C, 7 days at 2-8 °C and 1 year at -70 °C (freeze only once).

Preparation of samples

Dispense 1 mL of reagent R3 in a plastic or glass tube. Add 20 µL of well mixed blood sample (calibrator),

control or patient sample). Allow to stand for 5 minutes or until lysis is complete. Haemolysed samples are stable for 7 days at 2-8 °C.
On some analysers, hemolysis can be done on-board.

REFERENCE VALUES (4,5)

Recommendations for the follow-up of diabetic patients have been published in 2008 by ADA (American Diabetes Association).

	NGSP/DCCT (%)	IFCC (mmol/mol)
Non-diabetics	4.0 - 6.0	20 - 42
Goal under therapy (adults, except pregnant women)	< 7.0	< 53
Level 1	80	4.2
Level 2	80	7.1
Level 3	80	10.4

However, a level as close as possible to normal level (6.0 % or 42 mmol/mol) is recommended while avoiding hypoglycaemia. Less stringent goals may be applied in some cases (history of severe hypoglycaemia, children).

Note : The quoted range should serve as a guide only. It is recommended that each laboratory verifies this range or establishes a reference interval for the intended population.

PROCEDURE

For ELITech Clinical Systems Selectra Analyzers, applications are available on request.

Wavelength 660 nm

Temperature: 37 °C

	CALIBRATION	TEST
Reagent R1	225 µL	225 µL
Read absorbance (A1)		
Hemolysed calibrator *	6 µL	-
Hemolysed sample *	-	6 µL
Mix and after 4 minutes and 43 seconds, add :		
Working reagent R2	75 µL	75 µL
Mix and after 4 minute of incubation, measure the absorbance (A2)		

* Note : On some analysers, hemolysis can be done on-board.

With Selectra TouchPro software, use the application included in the barcode available at the end of this insert.

CALCULATION

Rate of HbA1c is calculated from a non linear calibration curve obtained from four calibrators of different levels and a zero point:

$$HbA1c = f(A2-A1)$$

CALIBRATION (5)

For calibration use ELITech HbA1c Calibrator Set, ref. HBAC-0043.

Add a zero point (NaCl, 9 g/L solution).

The values obtained with ELITech HbA1c reagent have been standardized from NGSP/DCCT (%) System. These values can be converted by calculation to traceable values from IFCC reference method (mmol/mol).

Following equations can be used to make this unit conversion :

$$HbA1c = f(A2-A1)$$

Calibrated Hemolysed: No significant interference up to 29.5 mg/dL (513 µmol/L).

Triglycerides: No significant interference up to 3046 mg/dL (34.42 mmol/L).

Acetyl salicylic acid: No significant interference up to 200 mg/dL.

Rhume factor: No significant interference up to 1000 IU/mL.

Ascorbic acid: No significant interference up to 20 mg/dL.

Interference due to Vitamin E has not been assessed.

- Specificity

Hemoglobin variants

Samples containing hemoglobin variants obtained from a NGSP reference lab were tested on ELITech Clinical Systems HbA1c. A no significant interference is defined by a recovery within +/-10% of NGSP reference value.

- HbC, HbS, HbD, HbE and HbA2: No significant interference.

- Samples containing HbF levels of 4% HbF up to 20%

HbF were tested. Samples up to 6% HbF showed no significant interference. A sample containing 7% HbF showed a bias of -12.9 %, a sample containing 20% HbF showed a bias of -21.5 %.

Chemically modified hemoglobin

The following HbA1c levels were tested: 6 % and 9 %. A no significant interference is defined by a recovery within +/-10% of the initial value.

- Carbamylated hemoglobin: No significant interference up to 10.0 mmol/L of added sodium cyanate.

- Acetylated hemoglobin: No significant interference up to 10.0 mmol/L of added aspirin.

- Labile HbA1c: No significant interference up to 1000 mg/dL of added glucose.

- It has been reported that results may be unreliable with patients who suffer from alcoholism (formation of acetaldehyde hemoglobin).⁽¹⁰⁾

- Many other substances and drugs may interfere. Some of them are listed in reviews published by Young.⁽¹³⁻¹⁴⁾

The exact range depends on the calibrator set used.

If the result is higher than 16 % (> 151 mmol/mol), report result as > 16 % (> 151 mmol/mol), do not dilute the sample.

Disposal of all waste material should be in accordance with local, state and federal regulatory requirements.

PERFORMANCE DATA at 37 °C on ELITECH Clinical Systems Selectra ProM Analyzers

- Measuring range

Determined according to EP6-A protocol of CLSI⁽⁶⁾, the reagent is linear from 2.5 % to 16 % (4 to 151 mmol/mol).

The exact range depends on the calibrator set used.

If the result is higher than 16 % (> 151 mmol/mol), report result as > 16 % (> 151 mmol/mol), do not dilute the sample.

Dispense 1 mL of reagent R3 in a plastic or glass tube. Add 20 µL of well mixed blood sample (calibrator),

- Limit of Blank (LoB) and Limit of Detection (LoD)
Determined according to CLSI⁽⁷⁾ EP17-A protocol.

LoB = 0.7 %

LoD = 0.8 %

- Precision

Determined according to EP5-A2 protocol of CLSI.⁽⁸⁾

	Mean	Within-run Total
n	%	mmol/L/mol
CV (%)		
Level 1	80	4.2
Level 2	80	7.1
Level 3	80	10.4

- Correlation

A comparative study has been performed according to CLSI EP9-A2 protocol⁽⁹⁾ between an ELITech Clinical Systems Selectra ProM Analyzer and another FDA-approved system equipment (equivalent method) on 40 human blood samples.

The sample concentrations were between 4.0 and 13.9 % (20 and 128 mmol/mol). The parameters of linear regression are as follows :

$$\text{Correlation coefficient: } r = 0.998$$

Linear regression: $y = 1.007x + 0.0\%$

- Limitations/Interferences^(1,2,10)

- This device should not be used for the diagnosis of diabetes mellitus and in judging day-to-day glucose control, and should not be used to replace daily home testing of blood glucose.

The following HbA1c levels are based on the Selectra ProM system. The reagent is linear from 2.5 % to 16 % (4 to 151 mmol/mol).

Interpretation of HbA1c values is based on erythrocytes lifespan; any change in this lifespan induces HbA1c rates which are not interpretable with reference values presented here (hemolytic diseases, recent significant blood loss, pregnancy).

- Not report results outside of the usable range.

Studies have been performed to determine the level of interference from different compounds according to CLSI EP7-A2 protocol⁽¹¹⁾ and SFBC recommendations⁽¹²⁾.

The following HbA1c levels were tested: 7 % and 10 %.

A no significant interference is defined by a recovery within +/-8% of the initial value.

Unconjugated Bilirubin: No significant interference up to 30 mg/dL (513 µmol/L).