## TRIGLYCERIDES MONO SL NEW

☛ **Références/Referencias/Referencias/Referências:**  
 TGML-0425  
 TGML-0515  
 TGML-0700  
 TGML-0427  
 TGML-0497  
 TGML-0517  
 TGML-0707

**Composition du coffret/ Kit composition/ Composición del kit/ Conteúdo da embalagem :**  
**R** 6 x 50 mL  
**R** 6 x 100 mL  
**R** 4 x 250 mL  
**R** 6 x 50 mL + **Std** 1 x 5 mL  
**R** 1 x 100 mL + **Std** 1 x 5 mL  
**R** 6 x 100 mL + **Std** 1 x 5 mL  
**R** 4 x 250 mL + **Std** 1 x 5 mL

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### BIBLIOGRAPHIE/BIBLIOGRAPHY BIBLIOGRAFIA/BIBLIOGRAFIA

- Naito, H.K., *Coronary Artery Disease and Disorders of Lipid Metabolism. Clinical Chemistry: Theory, Analysis, Correlation*, 4th Ed., Kaplan, L.A., Pesce, A.J., Kazmierczak, S.C. (Mosby, Inc. eds. St Louis USA), (2003), 603.
- Expert Panel on Detection, Evaluation, and Treatment of High Blood Cholesterol in Adults (*Adult Treatment Panel III*), *Executive Summary of the Third Report of the National Cholesterol Education Program (NCEP)*. *JAMA*, (2001), **285**, 2486.
- Fossati, P., Prencipe, L., *Serum triglycerides determined colorimetrically with an enzyme that produces hydrogen peroxide. Clin. Chem.*, (1982), **28**, 2077.
- Tietz, N.W., *Clinical guide to laboratory tests*, 3rd Ed. (W.B. Saunders eds. Philadelphia USA), (1995), 610.
- Evaluation of the Linearity of the Measurement of Quantitative Procedures: a Statistical Approach: Approved Guideline*. CLSI(NCCLS) document EP6-A (2003), **23** (16).
- Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline-Second Edition*. CLSI (NCCLS) document EP5-A2 (2004), **24** (25).
- Method Comparison and Bias Estimation Using Patient Samples; Approved Guideline-Second Edition*. CLSI (NCCLS) document EP9-A2 (2002), **22** (19).
- Interference Testing in Clinical Chemistry ; Approved Guideline-Second Edition*. CLSI (NCCLS) document EP7-A2 (2005), **25** (27).

9. Berth, M. & Delanghe, J. *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature. Acta Clin Belg.*, (2004), **59**, 263.  
 10. Young, D.S., *Effects of preanalytical variables on clinical laboratory tests*, 2nd edition, AACCPress (1997).  
 11. Young D.S., *Effects of drugs on clinical laboratory tests*, 4th edition, AACCPress (1995).

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- Only for ref. **TGML-0250/TGML-0455**, used with Selectra TouchPro software.  
 - **See § PROCEDURE:** Contamination risk

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

GPO = Glycérol-3-phosphate oxydase  
 DHA-P = Dihydroxyacétone phosphate  
 4-AAP = Amino-4-antipyrine

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



**PRINCIPLE** <sup>(3)</sup>  
 Enzymatic determination of triglycerides according to the following reactions :

GPO = Glycerol-3-phosphate oxidase  
 DHA-P = Dihydroxyacetone-Phosphate  
 4-AAP = Amino-4-antipyrine

### COMPOSITION

**Reagent: R**  
 Good's buffer, pH 7.00  
 p-Chlorophenol 2,7 mmol/L  
 ATP 3,15 mmol/L  
 Amino-4-antipyrine (4-AAP) 0,31 mmol/L  
 Lipoprotein lipase ≥ 2 000 U/L  
 Glycerol kinase ≥ 500 U/L  
 Glycerol-3-phosphate oxidase ≥ 4 000 U/L  
 Peroxidase ≥ 500 U/L  
 Sodium azide < 0,1 %  
 Also contains magnesium salts, FAD and detergents for optimal performance.  
**Standard: Std** (Ref : TGML-0427/497/0517/0707)  
 Glycerol (triglycerides equivalent) 200 mg/dL  
 2,26 mmol/L  
 Sodium azide < 0,1 %

### MATERIALS REQUIRED BUT NOT PROVIDED

- CALI-0550 ELICAL 2
- CONT-0060 ELITROL I
- CONT-0160 ELITROL II
- Normal saline solution (NaCl 9 g/L).
- General Laboratory equipment.
- Do not use materials that are not required as indicated above.

### WARNINGS AND PRECAUTIONS

- These *in vitro* diagnostic devices (reagent and standard) are for professional use only.
- The reagent R and standard contain sodium azide which may react with lead or copper plumbing to form potentially explosive metal azides. When disposing of these reagents always flush with copious amounts of water to prevent azide buildup.
- Take normal precautions and adhere to good laboratory practice.
- Use clean or single use laboratory equipment only to avoid contamination.
- The standard should be immediately and tightly capped to prevent contamination and evaporation.
- For more information, Safety Data Sheet (SDS) is available on request for professional user.

### STABILITIES

Store at 2-8 °C and protect from light. Do not freeze. Do not use after expiration dates indicated on the vial labels.  
**On board stability :**  
 The on-board stability is specific for each analyzer. (Refer to § PERFORMANCE DATA).

### PREPARATION

The reagent and standard are ready to use.

### PRODUCT DETERIORATION

- The product should be clear. Cloudiness would indicate deterioration.
- Do not use the product if there is visible evidence of biological, chemical or physical deterioration.
- Do not use the product if the damages of packaging might have an effect on the product performances (leakages, pierced vial).

### SAMPLES

- Serum or lithium heparinized plasma from fasting patients (≥ 12 hours).
- Do not use icteric or hemolyzed samples.
- Do not use other specimens.
- Warnings and precautions**
- According to Good Laboratory Practice, sampling should be performed prior to the administration of drugs. Sampling could lead to false results if performed during or immediately after the administration of some drugs.
- Collect the samples in tubes and stoppers free of glycerol.

- Separate from cells within 2 hours.
- Storage and stability**
- Samples are stable 5 to 7 days if stored at 2-8 °C, 3 months at -15 to -20 °C and several years at -70 °C. Avoid repeated freezing and thawing.

### REFERENCE VALUES

The NCEP (American National Cholesterol Education Program) has established the following classification for triglycerides levels according to the risk of developing coronary heart disease:  
**Risk Classification :**

	Level (mg/dL)	Level ( mmol/L)
Normal	< 150	1.69
Borderline high	150 - 199	1.69 - 2.25
High	200 - 499	2.26 - 5.64
Very High	≥ 500	5.65

**Note :** The quoted range should serve as a guide only. It is recommended that each laboratory verifies this range or establishes a reference interval for the intended population.

### PROCEDURE

**Manual Procedure**  
 Wavelength : 505 nm  
 Optical path : 1 cm  
 Sample/reagent ratio : 1:100  
 Temperature: 37 °C  
 Read against reagent blank.

	CALIBRATION	TEST
<b>Reagent R</b>	1 000 µL	1 000 µL
<b>Standard/Calibrator</b>	10 µL	-
<b>Sample</b>	-	10 µL

Mix and read the absorbances (A) after an incubation of 10 minutes.

### Automatic Procedure

These reagents may be used in several automatic analyzers. For ELITech Selectra Analyzers, validated applications are available on request. For Selectra TouchPro software, use the application included in the barcode available at the end of this insert.

- Triglycerides SL and Triglycerides Mono SL New reagents can be contaminated by Cholesterol HDL SL 2G reagent.  
 In order to avoid contamination on Selectra ProM and ProXL, program incompatibilities as follows:

Software	Menu	Parameter
TouchPro	Probe incompatibilities	Link / Cholesterol HDL SL 2G – Acid Solution
Other	Needle incompatibility	Cholesterol HDL SL 2G <<HCl

For other Selectra Pro instruments: repeat any absurd results after programming a needle wash.

- Lipase SL reagent is strongly contaminated by Triglycerides SL reagent.

In order to avoid cuvette contamination on Selectra Pro instruments, program the following incompatibilities:

Software	Menu	Parameter
TouchPro	Test incompatibilities	Link / Triglycerides SL – Acid Solution
Other	Cuvette incompatibility	Triglycerides SL <<HCl

In order to avoid needle contamination on Selectra Pro instruments, do not program Lipase SL and Triglycerides SL in the same run. Ensure the instrument goes back to "stand-by" status before launching a run containing Lipase SL.

### CALCULATION

$$\frac{A \text{ Sample}}{A \text{ Calibrator/ Standard}} \times n \quad n = \text{Calibrator/ standard concentration}$$
  
**Conversion factor :** mg/dL x 0,0113 = mmol/L  
 mg/dL x 0,01 = g/L

### CALIBRATION

For the references TGML-0427/0497/0517/0707: For calibration, use either multiparametric calibrator ELICAL 2 or Triglycerides Standard 200 mg/dL.  
 For the references TGML-0250/0455/0425/0515/0700: For calibration, use multiparametric calibrator ELICAL 2.

Concentration values of Triglycerides Standard 200 mg/dL and multiparametric calibrator ELICAL 2 are traceable to the reference method ID-GC-MS (Isotope Dilution - Gas Chromatography - Mass Spectrometry).

**Calibration frequency :** The calibration is specific for each analyzer. (Refer to § PERFORMANCE DATA).

### QUALITY CONTROL

To check the accuracy of assays, control sera such as ELITROL I and ELITROL II should be used. These controls must be performed and validated before the patient samples are assayed. The control frequency must be at least once a day, after each calibration and should be adapted to Quality Control procedures of each laboratory and the regulatory requirements. Results should be within the defined ranges. If values fall outside of the defined ranges, each laboratory should take corrective measures. Quality control materials should be used in accordance with local guidelines.

### WASTE MANAGEMENT

Disposal of all waste material should be in accordance with local, state and federal regulatory requirements.

### PERFORMANCE DATA at 37 °C on ELITech Clinical Systems Selectra ProM Analyzers

**Measuring range**  
 Determined according to CLSI<sup>®</sup> EP6-A protocol, the measuring range is from 30 to 1000 mg/dL (0.34 to 11,30 mmol/L).  
 Samples having greater concentrations should be diluted 1:5 with NaCl 9 g/L solution and re-assayed. This procedure extends the measuring range up to 5000 mg/dL (56.50 mmol/L).

For users with Selectra TouchPro software, the "dilute" function performs the sample dilution automatically. Results take the dilution into account.

### Precision

Determined according to CLSI EP5-A2 protocol<sup>®</sup>.

	Mean	Within-run	Total
<b>Low level</b>	80	44	0.50
<b>Medium level</b>	80	131	1.48
<b>High level</b>	80	267	3.02

### Correlation

A comparative study has been performed between an ELITech Clinical Systems Selectra ProM Analyzer and another FDA-approved system equipment (enzymatic and colorimetric method) on 99 human serum samples according to CLSI EP9-A2 protocol<sup>®</sup>. The values were between 30 and 957 mg/dL (0.34 and 10,81 mmol/L).  
 The parameters of the linear regressions are as follows :  
 Correlation coefficient: (r) = 0,999  
 Linear regression: y = 1,019 x + 1 mg/dL (0.10 mmol/L)

**Limitations/Interferences**  
 - Do not report results outside of the usable range.

- Studies have been performed to determine the level of interference from different compounds according CLSI EP7-A2 protocol<sup>®</sup>. Recovery is within ±10% of initial value of Triglycerides concentration of 133 mg/dL and 266 mg/dL.  
**Unconjugated Bilirubin:** No significant interference up to 15 mg/dL (257 µmol/L).  
**Conjugated Bilirubin:** No significant interference up to 5,9 mg/dL (101 µmol/L).  
**Hemoglobin :** No significant interference up to 125 mg/dL.  
**Uric acid :** No significant interference up to 24,2 mg/dL (1440 µmol/L).  
**Ascorbic acid :** No significant interference up to 2,0 mg/dL. Concentrations above the therapeutic levels will interfere and cause erroneous results.  
**Methyl-dopa :** No significant interference up to 1,0 mg/dL.

- In very rare cases, monoclonal gammopathies (multiple myeloma), in particular IgM type (Waldenström's macroglobulinemia) can cause unreliable results.<sup>(9)</sup>  
 - Results can be falsely lowered by significant levels in sample of NAC (N-Acetyl-Cysteine), NAPQI (metabolite of acetaminophene (paracetamol)) or metamazole.  
 - Many other substances and drugs may interfere. Some of them are listed in reviews published by Young.<sup>(10-11)</sup>

- The results of this assay should only be interpreted in conjunction with other diagnostic test results, clinical findings and the patient's medical history.

**On board stability/Calibration frequency**  
**On Board Stability:** 28 days  
**Calibration frequency:** 14 days  
 Recalibrate when reagent lots change, when quality control results fall outside the established range, and after a maintenance operation.

*These performances have been obtained using ELITech Selectra ProM analyzer. Results may vary if a different instrument or a manual procedure is used. The performances of applications not validated by ELITech are not warranted and must be defined by the user.*

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